Proline®

Instruction Manual

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
Mode d'emploi


Manual de Instrucciones

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Instruction Manual.....	1
Bedienungsanleitung.....	23
Mode d'emploi.....	45
Manual de Instrucciones.....	67
Istruzioni d'impiego.....	89

(*) PP-machinery, complete, includes: DRV-board, motor, membrane and linear actuator,
 (**) CPU-assembly, complete, includes: CPU-board, LCD, Keyboard + connectors.

cat.no.: 712726
 cat.no.: 712759

CONTENTS

1. YOUR NEW BIOHIT PROLINE XL	2
2. UNPACKING & PREPARING THE PIPETTOR FOR USE	2
2.1. Installing the Charging Base.	3
2.2. Charging the XL	3
2.3. Battery and the LO BAT function	4
2.4. Electrical Specifications	5
3. XL MATERIALS	6
4. XL DESCRIPTION	7
4.1. Controls.	7
4.2. Start Buttons	7
4.3. Display.	7
4.4. Direction Symbols	8
4.5. Nose cones and adapters.	8
4.6. Filters.	9
5. PROGRAMMING THE XL	9
5.1. Mode Selection and Mode Recall.	9
5.2. Setting Speeds	9
5.3. Setting the Pipetting range (the pipette/tip size to be used)	10
5.4. Pipetting Mode (P)	11
5.5. Multi-Dispensing Mode (d)..	11
5.6. Sequenced Dispensing Mode (Sd)... ..	12
5.7. Continuous Pipetting Mode (C)	13
6. PIPETTING RECOMMENDATIONS	13
6.1. Dispensing with Blow-out.	13
6.2. Dispensing without Blow-out... ..	14
6.3. Other Recommendations	14
7. STORAGE	14
8. CALIBRATION	15
8.1. Testing the Performance	16
9. MAINTENANCE	17
9.1. Decontamination.	17
9.2. Changing the filter, silicon adapter or nose cone	17
10. TROUBLE-SHOOTING	17
11. WARRANTY INFORMATION	19
12. SPECIFICATIONS	19

1. YOUR NEW BIOHIT PROLINE XL

Your new Biohit Proline XL, an electronic pipettor/pipetting controller, has been designed for maximum performance and accuracy with ergonomic and user-friendly operation. The microprocessor-based system reduces the possibility of human error and instrument contamination by controlling all piston movements. Its light-weight and ergonomic controls take the effort out of pipetting and help reduce the risk of repetitive strain injuries (RSIs) that are common in manual pipetting. The swivel head of the XL also rotates to an optimal position thus reducing strain to the wrist during use.

The Biohit Proline XL operates on the air displacement principle and is equipped with three different nose cones and adapters to accommodate all standard graduated or non-graduated glass or plastic pipettes as well as Biohit Proline 5 ml (780200 and 780300) and 10 ml (780310) polypropylene tips.

2. UNPACKING & PREPARING THE PIPETTOR FOR USE

The Biohit Proline XL package contains the following items:

- ☀ XL Pipettor/Pipetting aid with 25 ml nose cone (with silicon adapter)
- ☀ 10 ml nose cone (with silicon adapter) and 10 ml tip
- ☀ 5 ml nose cone (with silicon adapter) and 5 ml tip
- ☀ Charging base
- ☀ AC-adaptor
- ☀ Filters
- ☀ Instructions for use
- ☀ Performance certificate in accordance with DIN E-12650-7

Make sure that all items are included and that no damage has occurred during shipment.

2.1. Installing the Charging Base

1. Connect the AC-adaptor into the connector at the rear of the charging base (Fig. 1).
2. Before connecting the AC-adaptor to the electrical outlet make sure that the voltage setting corresponds to the voltage from your electrical outlet and that the AC-adaptor is compatible with your electrical outlet.
3. The green light in the front of the charging base indicates it is in operation.

Fig. 1

Note: The charging base should only be used with the AC-adaptor provided. Do not use electrical outlets that are timed or switched.

2.2. Charging the XL

The charging base operates via the spring loaded metal contacts on the bottom of the XL charging base. The red ON/OFF switch is located at the rear of the XL (Fig. 2).

1. Switch ON the XL (the far right position) (Fig. 2).
2. Place the pipettor in the charging base. To ensure proper charging make sure that the contacts in the charging base and the bottom of the XL are always clean (Fig. 3).

Fig. 2

Fig. 3

3. If the pipettor is new or the battery is low, keep the XL in the base for 14 hours to fully charge it before continuing use.
 4. An "r" sign will appear on the display when the XL is ready for use. Once you press either of the START buttons the XL is ready for both default pipetting at maximum volume and for program changes.
- ☀ The XL will admit a charge with the switch in either ON or OFF positions. However, if the XL is charged with the switch in the OFF position, the battery will not obtain a maximum charge as only a low current charge is applied.
 - ☀ Leave the pipettor ON in the charging base to retain user-selected settings as switching it OFF will reset the unit to default settings.
 - ☀ If the XL is left ON uncharged for several days, the display will go blank and there will be no response from the keyboard or START buttons, as the battery voltage has fallen below the required operating level.

2.3. Battery and the LO BAT function

The XL has a software/hardware-based battery capacity control system. During the recharging process five rectangular symbols will scroll on the display. When the battery is completely charged, a *FULL* signal will appear on the display and the XL will be ready for at least 500 dispensing cycles! The XL can be used without a *FULL* recharge, but we suggest that you allow the unit to reach a *FULL* charge whenever possible.

As the battery voltage decreases during use, the amount of rectangles on the display will decrease respectively. When the battery is close to empty, the five rectangles on the display will start to blink, indicating the need to recharge the unit (Fig 4). When the blinking starts, the unit can still be used for a short period of time. However, the XL should be returned to the base for recharging as soon as possible.

Note: Always keep the power switch in the ON position as the XL will then automatically switch into a "sleeping mode". This mode guarantees minimum energy consumption when the XL is not in operation.

Fig. 4

2.4. Electrical Specifications

Battery


- ☀ Rechargeable NiCd battery
- ☀ Charging time max 14 hours for a completely empty battery
- ☀ The battery pack has a built in short-circuit protection to aid in the prevention of short circuits

AC-Adaptor

- ☀ Input voltage and main plug according to local requirements
- ☀ Output voltage 9 VDC

Note: Do not use any other than the original manufacturer's AC-adaptor

3. XL MATERIALS

Note: The chemical resistance of the cylinder material is limited. For the demanding solvent applications, a more durable HDPE cylinder is available as a spare part (cat. no. 712757).

4. XL DESCRIPTION

The control and programming of the Biohit Proline XL are done using the keyboard and display shown in detail in figure 5.

4.1. Controls

Fig. 5

- S** Changes speeds
- M** Moves to next mode
- ▲** Increases volume or speed
- ▼** Decreases volume or speed
- E** Confirms setting changes
- *** Displays range and confirms setting changes

4.2. Start Buttons

There are two START buttons located at the front of the handle. In the continuous mode (C), the upper button enables aspiration for as long as the button is held down (within the operating range). The lower button allows for dispensing. Furthermore, the aspiration/dispensing speed can be regulated by applying a different amount of finger pressure to the START buttons. When the programmable modes (P, d, Sd) are employed, either of the START buttons can be used to activate the unit once all the required information has been programmed. Only a light touch is needed to activate a function on these modes. If the programmed function must be stopped, e.g. when the programmed volume exceeds the volume of the pipette/tip in use, pressing both buttons simultaneously will interrupt the operation and a **STOP** signal will appear on the display. Pressing either of the buttons once more will reset the unit, i.e. empty the pipette/tip.

4.3. Display

The left display is the status indicator. It informs the user about parameters to be programmed, functions to be performed, the number of dispensings available, and the remaining capacity of the battery.

The right display is used for programming and displaying the various volumes needed in different operating modes. It also contains the piston direction symbols.

4.4. Direction Symbols

These symbols indicate the direction in which the piston moves upon pressing the START button. The small UP ARROW on the display means that the next function is aspiration. The small DOWN ARROW in turn indicates dispensing. Both functions will be carried out in accordance with the selected operating mode.

4.5. Nose cones and adapters

The XL is provided with three different pipette/tip nose cones/adapters. The 25 ml nose cone/adapter allows the use of all standard plastic and glass pipettes (graduated and non-graduated) including Pasteur pipettes. The 5 and 10 ml nose cone/adapters are intended to be used with Biohit Proline 5 and 10 ml tips, accordingly.

To change the nose cones, filters or silicon adapters, first press the nose cone lock release button (Fig. 6), rotate the nose cone counter-clockwise and pull out the nose cone and the silicon adapter. Insert the adapter/nose cone of your choice and twist the nose cone gently clockwise. The silicon adapters enable the pipettes/tips to fit securely in place, thus creating proper vacuum levels and reducing leakage. Both the nose cones and silicon adapters are autoclavable (121°C, 1 atm).

Fig. 6

Note: It is recommended to change the silicon adapter every six months, or when necessary, to prevent a loss of vacuum. Spare adapters can be ordered from the local Biohit distributor.

Cat. No	Product	Qty.
712893	25 ml adapter for XL	1 pc
712911	10 ml adapter for XL	1 pc
712910	5 ml adapter for XL	1 pc

4.6. Filters

Removing the nose cone and adapter exposes a hydrophobic filter (0.45 micron) located at the top of the nose piece. The filter protects the XL from accidental over-aspirations and we strongly recommend to use only the original Biohit filters. Both sterile and non-sterile filters are available.

Cat. No	Product	Qty.
712912	Filter for Biohit Proline XL, non-sterile	5 pcs
712913	Filter for Biohit Proline XL, sterile	1 pcs

Note: The XL can not be operated without a filter. Wet filters, due to overfilling will slow or stop fluid movement. If overfilling occurs, change the filter immediately by simply pulling it out and placing a new filter, with the blue side upwards, into the nose piece.

5. PROGRAMMING THE XL

Programming is done using the six-button keyboard and the LC-display. There are four different operating modes, offering special functions and varying speeds for your selection.

5.1. Mode Selection and Mode Recall

1. Press **M** repeatedly to view the modes available.
2. Press **E** when the desired mode is displayed. The XL is ready for operation on the selected mode.

Note: A mode can only be activated when the piston is in its home position (arrow up sign is lit).

5.2. Setting Speeds

1. Press **S** to display the current aspiration speed.
2. Press **▲** or **▼** until the desired aspiration speed is displayed (5 = Fast and 1 = Slow).
3. Press **E** to confirm the speed selection. The display shows the current dispensing speed.
4. Press **▲** or **▼** until desired dispensing speed is displayed (5 = Fast and 1 = Slow).
5. Press **E** to confirm the speed selection.

The default speed for both aspiration and dispensing for each pipetting range is as follows:

Pipetting range	Default speed	Pipetting range	Default speed
25 ml	5	2 ml	2
10 ml	4	1 ml	1
5 ml	3		

Note: The speed can not be changed in the middle of the aspiration or dispensing cycles.

5.3. Setting the Pipetting range (the pipette/tip size to be used)

To avoid accidental over aspiration the XL can be programmed to perform on a selected range when using smaller than 25 ml pipettes/tips. There are designated limits for 1, 2, 5, 10 and 25 ml tips. The selected limit remains active for all the operating modes, e.g. in the C-mode the piston movement stops if the aspiration goes beyond the preset volume limit, and in the P-, d-, and Sd-modes the programming can be done only within the allowed range.

1. Press ***** to display the current range.
2. Press **▲** or **▼** until the desired range is displayed.
3. Press **E** to confirm range selection.

Note: The chosen volume range remains active for all the operating modes unless readjusted.

Maximum volume	Range (P-mode)	Range (d-mode)	Range (Sd-mode)	Range (C-mode)
1 ml	0.1-1.0 ml	0.1-0.5 ml	0.1-1.0 ml	0-1.0 ml
2 ml	0.2-2.0 ml	0.1-1.0 ml	0.1-2.0 ml	0-2.0 ml
5 ml	0.5-5.0 ml	0.1-2.5 ml	0.1-5.0 ml	0-5.0 ml
10 ml	1.0-10.0 ml	0.5-5.0 ml	0.5-10.0 ml	0-10.0 ml
25 ml	2.0-25.0 ml	1.0-12.5 ml	1.0-25.0 ml	0-25.0 ml

5.4. Pipetting Mode (P)

The pipettor performs ordinary pipetting with blow-out.

1. Press **M** to display .
2. Press **E** to confirm the mode change.
3. Select the desired pipetting volume by pressing  to increase, and  to decrease the volume.

Note: When holding down the  or the  button, the volume display starts to scroll faster.

4. Press **E** to confirm the selection.
5. Position the tip to aspirate and press either of the START buttons.
6. Position the tip to dispense. Press either of the START buttons. The tip is emptied with a blow-out and is ready for the next pipetting.

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The XL performs repetitive dispensings of a selected volume. During this operation, the desired volume plus an automatically selected excess volume is aspirated into the tip. The excess volume is needed to ensure equal operating conditions for each dispensing.

1. Press **M** to display .
2. Press **E** to confirm the mode change.
3. The volume and number of dispensing aliquots appear in the display. If this mode has been used before the previous settings appear in the display. Press **E** to start the operation or ***** to change the volume.
4. Press  or  until the desired multi-dispensing volume is displayed.
5. Press ***** to confirm the selection.
6. Press  or  until the desired number of aliquots is displayed.
7. Press **E** to confirm the selection.
8. Position the tip to aspirate and press either of the START buttons.
9. The * sign and DOWN ARROW are lit to indicate the reset function. Position the tip to discard priming excess and press either of the START buttons.

10. To dispense, position the tip, press either of the START buttons and repeat until the cycle is complete.

Note: To avoid dripping, the XL will automatically aspirate a small amount of air after each dispensing.

11. Finally, position the tip to discard any remaining excess and press the START button twice.

Note: The dispensing can be interrupted by pressing **E**. An ARROW DOWN and E appear on the display. Pressing either of the START buttons twice will empty the pipette/tip.

5.6. Sequenced Dispensing Mode (Sd)

A series of different volumes can be delivered in any desired order.

1. Press **M** to display **Sd**.
2. Confirm the mode change by pressing **E**.
3. The first dispensing volume appears in the display. If this mode has been used before the previous settings appears in the display. Press **E** to start the operation or ***** to change the volume.
4. Select the first desired dispensing volume onto the display by using the **▲** and **▼**.
5. Press ***** to confirm the desired dispensing volume.
6. To select the following dispensing volumes (up to xx) use the **▲** and **▼** buttons, always remembering to confirm each selection by pressing the ***** button.
7. To confirm your final selection press **E**.
8. Position the tip to aspirate and press either of the START buttons.
9. Position the tip to discard priming excess and press either of the START buttons again.
10. Position the tip to dispense and press either of the START buttons. Repeat the action until the cycle is complete.

Note: To avoid dripping, the XL will automatically aspirate a small amount of air after each dispensing.

11. Position the tip to discard any remaining excess and press START twice.

5.7. Continuous Pipetting Mode (C)

The aspiration and dispensing are controlled manually. The aspiration/dispensing continues as long as the corresponding button is held down (within the operating volume range). The C-mode is also provided with the blow-out function.

1. Press **M** to display .
2. Press **E** to confirm the mode change.
3. The home position, i.e. 0.00 ml, will be shown on the display.
4. Press the upper START button to begin aspiration. Keeping the button down sustains the function until the maximum volume of the chosen range has been reached. The function can be paused or restarted at any moment by manipulating the START button.
5. Press the lower START button to begin dispensing. Keeping the button down sustains the function until the home position has been reached.
6. If the blow-out function is needed to empty the pipette/tip completely, a double click of the START button will initiate the blow-out function and the piston is returned back to the home position.

- ☀ The harder the buttons are pressed the faster the aspiration/dispensing.
- ☀ While the piston is moving, a C and a dashed line (C - - - -) appears on the display.
- ☀ If the maximum volume is aspirated, the XL will automatically stop and the aspirated volume will appear on the display.

6. PIPETTING RECOMMENDATIONS

Several different liquid handling procedures are available when the different operating modes and various special functions are utilised. Two of the modes (P and C) feature an automatic blow-out and others (d and Sd) leave the excess liquid in the tip. The recommendations below should be followed to ensure optimal performance.

6.1. Dispensing with Blow-out

The pipetting mode (P), has an automatic blow-out function, and the continuous mode (C), a manual blow-out. In both modes the blow-out function is followed by an automatic return of the piston to the "home" position. To avoid accidental aspiration of the liquid back into the tip, it is recommended that dispensing is always conducted above the liquid surface.

6.2. Dispensing without Blow-Out

The blow-out function cannot be employed on the dispensing modes (d, Sd). Therefore, it is recommended that dispensing is always performed with the pipette/tip set against the wall of the container.

Note: To avoid dripping, the XL will automatically aspirate a small amount of air after each dispensing.

6.3. Other Recommendations

- ✱ Hold the XL vertically when aspirating the liquid and place the tip only a few millimetres into the liquid.
- ✱ When using automatic modes, pre-rinse the pipette/tip by filling and emptying the tip for five times before aspirating the actual sample. This is important especially when dispensing liquids which have a viscosity and density different from water.
- ✱ When pipetting liquids at a temperature other than ambient, pre-rinse the pipette/ tip several times before use.
- ✱ Do not rest the pipettor on its side, as this could cause contamination of the nose cone and filter. Always place the XL in the charging base when it is not in use.
- ✱ Do not drop the XL, the charging base or the AC-adaptor as this may damage these devices.
- ✱ Avoid exposing the unit to extreme temperature changes, humidity and dust.
- ✱ Avoid rough handling, moderate pressure is all that is required when using the keyboard.

7. STORAGE

It is recommended that the XL is kept in its charging base with the switch in the ON position when it is not in use.

8. CALIBRATION

The calibration of each Biohit Proline XL has been factory checked and certified at 22°C using distilled water according to DIN 12650-7. Specifications of the XL are guaranteed only when it is used with genuine Biohit Proline tips. Due to the varying nature of different manufacturer's serological pipettes, no specifications for serological pipettes, others than 25 ml are given. The ranges of variation of serological pipettes are listed in the Specifications section (12).

The Biohit Proline XL can be re-calibrated in P-mode on volume ranges 5, 10, and 25 ml. Re-calibration cannot be performed on volume ranges 1 and 2 ml. The re-calibration can be performed in one selected volume at a time in 0.05 ml increments as follows.

1. Press **M** to display **P**.
2. Confirm the mode change by pressing **E**.
3. Select the desired pipetting volume to be calibrated by using **▲** to increase and **▼** to decrease.
4. Press **E** to confirm the selection.
5. Press **E** and **▲** simultaneously to display the volume to be calibrated.

Note: When pressing both buttons the display will be blank except for the "ml" sign. Releasing the buttons will display the volume to be calibrated.

6. Use **▲** to increase or **▼** to decrease the steps (in 0.05 ml steps) for the selected volume (e.g. for 25.00 ml volume in display program the pipettor to give out 25.05 ml or 24.55 ml).
7. Press **E** to confirm the selection.
8. Display will show **PC** and the original volume (which is now calibrated to give out the new volume).
9. Start pipetting.

Note: Changing the volume or mode will reset the calibration to default settings!

8.1. Testing the Performance

1. Carefully fit the pipette/ tip onto the nose cone.
2. Pre-rinse the tip with distilled water by pipetting the selected volume five times.
3. Carefully aspirate the liquid, keeping the pipettor vertical.
4. Pipette distilled water into a tared container and record the weight registered in mgs. Repeat at least ten times. Use an analytical balance with a readability of 0.01 mgs.
5. Compare the results with the corresponding specifications (section 12). If the mean value of 10 readings falls within the indicated range, the XL is ready for use.

- ☀ Weighing should take place at 20-25°C, constant to 0.5°C.
- ☀ Avoid drafts.
- ☀ The distilled water, weighing vessel, XL and pipettes/ tips must be at the same temperature.
- ☀ To calculate the volume, divide the weight of the water by its density (at 20°C: 0.9982). This method is based on DIN 12650-7.

9. MAINTENANCE

The Biohit Proline XL requires regular cleaning to ensure trouble-free operation. Use a soft lint free cloth lightly moistened with a mild detergent to clean the outer surface of the pipettor. DO NOT AUTOCLAVE. Only the nose cones and the silicon adapters are autoclavable (121°C, 1 atm).

Note: The pipettor must be turned off prior to servicing!

9.1. Decontamination

To decontaminate the XL, lightly spray the outside of the XL with the Biohit Proline Biocontrol decontamination spray (cat.no 724003, 500 ml or 724004, 5 l) and wipe dry with a clean lint free cloth. The nose cones and the silicon adapters (see below) can also be placed into a beaker containing Biohit Proline Biocontrol and left there for 30 minutes for complete decontamination. The parts should then be rinsed with distilled water and dried.

9.2. Changing the filter, silicon adapter or nose cone

To change nose cones, filters or silicon adapters, first press the nose cone lock release button (4.5 Nose cones and adapters, Fig. 6), rotate the nose cone counter-clockwise and pull out the nose cone, the silicon adapter and the filter when necessary. Insert the new filter with the blue side upwards and the adapter/nose cone of your choice and twist the nose cone gently clockwise. Both the nose cones and silicon adapters are autoclavable (121°C, 1 atm).

10. TROUBLE-SHOOTING

The Biohit Proline XL has a built-in monitoring program to control the precision and performance of each pipetting action. If the error message "Er" appears on the display, this means the XL has been unable to perform the attempted action properly. If you should receive the error message please do the following:

(As this procedure will empty the tip, it is recommended that you remove the tip before resetting the XL.)

1. Clear the error message from the display by pressing **E**.
2. Place the XL in the charging base for 15 minutes.
3. Press START, which will return the piston to its home position.

- ☀ Occasional "Er" situations can occur in the case that main outlets have been switched off. They can also take place if the operating switch has been in the OFF position during charging.
- ☀ Repeated occurrence of the the "Er" message is caused by internal error, which is due to a failure in the execution of the stroke path. In this case the pipettor will need to be returned to the local Biohit distributor for repair.

Problem	Possible cause	Solution
Droplets left inside	Unsuitable pipette/tip	Use original Biohit tips or pipette/tips specified by the manufacturer
Leakage or pipetted volume too small	Pipette/tip incorrectly attached	Reattach the pipette/tip
	Unsuitable pipette/tip	Use original Biohit tips or pipette/tips specified by the manufacturer
	Foreign particles between adapter and pipette/tip	Clean the adapter
	Silicon adapter worn out	Change the adapter
XL out of given specs	Check the pipette/tip in use	Change the pipette/tip
	Instrument damaged	Send for service
XL blocked, aspirated volume too small	Filter contaminated	Change the filter
	Adapter contaminated	Clean or change the adapter
Continuous error message	Instrument damaged	Send for service

11. WARRANTY INFORMATION

The Biohit Proline XL is warranted for one year (the battery is warranted for three months) against defects in materials and workmanship. Should it fail to function at any time, please contact your local representative immediately. The warranty will not cover defects caused by normal wear or by handling or using the pipettor in a manner contrary to the instructions given in the manual.

Each Biohit Proline XL is tested before shipping by the manufacturer. The Quality Assurance Procedure is your guarantee that the Biohit Proline XL you have purchased is ready for use.

Each Biohit Proline XL is CE-marked, fulfilling the requirements of the EMC standards EN 55014, 1993 and EN 55104, 1995.

12. SPECIFICATIONS

The Biohit Proline XL has been quality tested according to DIN 12650. The performance control according to DIN 12650 involves gravimetric testing of each XL with distilled water (quality 3, DIN ISO 3696) at 22°C using the manufacturer's original tips or pipettes specified by the manufacturer.

Note: The manufacturer's specifications below should be used as guidelines when establishing your own performance specifications in accordance with DIN 12650.

The performance for Biohit Proline 10 ml and 5 ml polypropylene tips:

Pipette size	Mode	Test volume	Inaccuracy	Imprecision
10 ml	P-mode	10 ml	0.50 %	0.15 %
10 ml	P-mode	1 ml	2.40 %	0.30 %
10 ml	d-mode	10 x 1 ml	1.00 %	1.30 %
5 ml	P-mode	5 ml	1.00 %	0.20 %
5 ml	P-mode	0.5 ml	4.20 %	0.60 %
5 ml	d-mode	10 x 0.5 ml	1.70 %	1.50 %

Due to the varying nature of different manufacturer's serological pipettes, no specifications for serological pipettes, others than 25 ml are given. The performance for 25 ml Sterilin serological pipette (factory test):

Pipette size	Mode	Test volume	Inaccuracy	Imprecision
25 ml	P-mode	25 ml	0.40 %	0.30 %
25 ml	P-mode	2.0 ml	1.40 %	0.60 %
25 ml	d-mode	10 x 2.5 ml	1.00 %	1.50 %

The performance of different manufacturer's serological pipettes vary due to the following properties:

- ☀ Length of the pipette (air space between the piston and liquid)
- ☀ Shape of the pipette (cone)
- ☀ Manufacturing technique
- ☀ Material (glass/plastic, quality of the plastic)

With most serological pipettes the precision is good but there may be big differences in accuracy. The following specifications are the guidelines for serological pipettes:

Pipette size	Mode	Test volume	Inaccuracy	Imprecision
25 ml	P-mode	25 ml	0.40 %	0.30 %
25 ml	P-mode	2.5 ml	2.40 %	1.00 %
25 ml	d-mode	10 x 2.5 ml	2.00 %	1.50 %
10 ml	P-mode	10 ml	0.80 %	0.30 %
10 ml	P-mode	1 ml	2.50 %	2.30 %
10 ml	d-mode	10 x 1 ml	2.50 %	2.50 %
5 ml	P-mode	5 ml	1.50 %	0.40 %
5 ml	P-mode	0.5 ml	4.50 %	2.30 %
5 ml	d-mode	10 x 0.5 ml	2.50 %	2.50 %
2 ml	P-mode	2 ml	2.60 %	1.50 %
2 ml	P-mode	0.2 ml	8.00 %	1.50 %
2 ml	d-mode	10 x 0.2 ml	3.00 %	3.00 %
1 ml	P-mode	1 ml	2.60 %	1.50 %
1 ml	P-mode	0.1 ml	8.00 %	1.50 %
1 ml	d-mode	10 x 0.1 ml	5.50 %	4.00 %

Biohit	Ex
Proline	20°C

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Due to Biohit's continuing R&D effort, specifications may change without prior notice. Biohit Proline innovations are covered by multiple patents including U.S. Pat. 5,343,769, U.S. Pat. 5,505,097, JP Pat. 3215860, JP Pat. 3240374 and patents pending. Biohit Proline® is a registered trademark of Biohit Plc.
Biohit S.A.
Biohit Deutschland GmbH
Biohit Japan Co., Ltd
Finnbio Ltd.
Biohit Ltd.