


VWR Signature™ Ergonomic High-Performance Pipettor



ENGLISH	1 – 8
DEUTSCH	9 – 18
FRANÇAIS	19 – 28
ESPAÑOL	29 – 38
ITALIANO	39 – 48

VWR Signature™ Ergonomic High-Performance Pipettor

VWR Signature™ Ergonomic High-Performance Pipettor

VWR Signature™ Ergonomic High-Performance Pipettor

ENGLISH

CONTENTS

- 1 - INTRODUCTION
- 2 - VOLUME SETTING
- 3 - METHOD OF PIPETTING
- 4 - RECOMMENDATIONS
- 5 - RECALIBRATION
- 6 - CLEANING AND STERILIZATION
- 7 - PIPETTOR KIT AND ACCESSORIES
- 8 - SPARE PARTS

1 - INTRODUCTION

The **VWR Signature™ Ergonomic High-Performance Pipettors** have been designed for use with microplates, loading electrophoresis gels and for other multisample applications. The pipettors enable precise and simultaneous delivery of 8 or 12 pre-set-volume doses of liquid.

VWR Signature™ Pipettors are produced in four ranges of volumes: from 1 µl to 10 µl, 5 µl to 50 µl, 20 µl to 200 µl and 50 µl to 300 µl.

The **VWR Signature™ Pipettors** have a digital counter displaying the set volume. The set volume is visible in the window situated in the hand grip of the pipettor. Volume is adjusted by turning the pipetting pushbutton (Fig. 1A2) or by turning the adjustment knob [Fig. 1B). The pipette design allows the user to lock the volume setting by pushing the locking ring (Fig. 1P) upwards. The position of the ring is indicated by the symbols located on the handle. The volume range of the pipettor is shown on the pushbutton (Fig. 1A1).

VWR Signature™ Pipettors should be used with polypropylene disposable tips, to ensure safety and accuracy (Fig. 1C). In order to protect the user against inadvertent contact with used tips, the pipettors are equipped with tip ejectors (Fig. 1D). The Pipette construction allows to rotate the multichannel module and to eject the tips through 360°, therefore, the most convenient position may be selected to fill the microplates. It is recommended to turn the module clockwise.

The **VWR Signature™ Pipettors** is a high quality instrument which offers excellent accuracy and precision. The

specifications for accuracy and precision given in the following table were obtained using **VWR** tips. These figures are only guaranteed when **VWR** tips are used.

8 channel

SPECIFICATIONS				
Model	Volume [µl]	Accuracy [%]	Precision [%]	Fit to tips µl
VWR ME10	Min 1	±8.0	≤ 6.0	10
	5	±4.0	≤ 2.0	
	Max 10	±2.0	≤ 1.2	
VWR ME50	Min 5	±4.0	≤ 2.5	200
	25	±3.0	≤ 1.2	
	Max 50	±1.6	≤ 0.6	
VWR ME200	Min 20	±3.0	≤ 1.5	200
	100	±1.5	≤ 0.8	
	Max 200	±1.0	≤ 0.6	
VWR ME300	Min 50	±1.6	≤ 1.5	300
	150	±1.2	≤ 1.0	
	Max 300	±1.0	≤ 0.6	

12 channel

SPECIFICATIONS				
Model	Volume [µl]	Accuracy [%]	Precision [%]	Fit to tips µl
VWR MT10	Min 1	±8.0	≤ 6.0	10
	5	±4.0	≤ 2.0	
	Max 10	±2.0	≤ 1.2	
VWR MT50	Min 5	±4.0	≤ 2.5	200
	25	±3.0	≤ 1.2	
	Max 50	±1.6	≤ 0.6	
VWR MT200	Min 20	±3.0	≤ 1.5	200
	100	±1.5	≤ 0.8	
	Max 200	±1.0	≤ 0.6	
VWR MT300	Min 50	±1.6	≤ 1.5	300
	150	±1.2	≤ 1.0	
	Max 300	±1.0	≤ 0.6	

These specifications are obtained in forward mode, using a gravimetric method with the temperature of the distilled water, tips and all other conditions stabilized between 19°C and 21°C. Number of measurements minimum 10. The values given include all components of error due to both normal handwarming and changing of the tips.

Performance tests: The **VWR Signature™ Pipettors** is calibrated based on EN ISO 8655. Performance can be verified by checking the pipette using the procedures outlined in the EN ISO 8655 standard.

The pipettor design enables the user to perform the recalibration process according to the rules presented in section 5.

VWR Signature™ Ergonomic High-Performance Pipettor

ENGLISH

2 - SETTING THE VOLUME

Volume is set using the pipetting pushbutton (Fig. 1A2) or adjustment knob (Fig. 1B), while the value is shown on the digital counter. The volume shown by the counter is represented by three digits which should be read from top to bottom. Typical readings are shown in the following table:

Model	Counter readings	Set volume	Basic degree
VWR ME10 VWR MT10	0 3 5	3.5 μ l	0.02 μ l
VWR ME50 VWR MT50	0 6 5	6.5 μ l	0.10 μ l
VWR ME200 VWR MT200	0 8 5	85.0 μ l	0.20 μ l
VWR ME300 VWR MT300	2 5 0	250.0 μ l	1.0 μ l

The volume adjustment can be performed when the locking ring is set in the lower position (Fig. 5A). When the desired volume is selected, the locking ring should be set in the upper position (Fig. 5B).

To attain the maximum accuracy, set volume must be approached from a higher value. If the desired value is lower than the previous value, simply use the knob in the pushbutton (Fig. 1A2) or the adjustment knob (Fig. 1B) to adjust the value down to 1/3 turn above the required setting. If the desired value is higher than the previous value, use the knob in the pushbutton or the adjustment knob to adjust the value to 1/3 turn above the required value. The setting should then be performed in the following manner:

With the volume set to 1/3 above the required setting, slowly turn the pipetting pushbutton knob or the adjustment knob to the required volume.

If the knob is accidentally turned too far, ie is set lower than the required value, the process must be repeated. The approach to the set volume must always be made in the order of decreasing value.

Following volume adjustment, set the locking ring into the upper position, thus locking the knob and preventing accidental volume change.

3 - METHOD OF PIPETTING

Fit the tips onto the cones of the multichannel manifold. While holding the pipettor in a vertical position, tips should be immersed in the liquid to a depth of 2-4 mm, and flushed once by drawing a dose of liquid and dispensing it out with slow and steady movement.

Next, still holding the pipettor vertically, press the pipetting button until the first resistance point is felt (Fig. 3A), and immerse the tips in the liquid to the depth of 2-4 mm. Releasing the pipetting button with a slow and uniform movement during 2-3 seconds, the liquid should be drawn into the tips (Fig. 3B) which should then be lifted above the liquid's surface. The pipettor should next be positioned at an angle of 10-45° in relation to inner walls of destination vessels and the tips should be emptied by pressing the pipetting button slowly until the first resistance point is felt (Fig. 3C). After waiting a second, the pipetting button should be pressed to the second resistance point in order to expel remaining liquid (Fig. 3D). Then the tips should be lifted out from the vessels while maintaining contact between the ends of the tips and the inner walls of the vessel until the pipetting button has been released. Finally, press the ejector button to separate the tips from multichannel manifold's cones and dispose of properly (Fig. 3F).

4 - RECOMMENDATIONS

The recommendations below will ensure maximum accuracy and precision from your **VWR Signature™ Pipettor**.

- When pipetting, the **VWR Signature™ Pipettor** should have the volume setting locked with the locking ring in the upper position.
- Do not draw liquids without tips fitted on the tip cones.
- Do not lay down pipettor with tips filled on the pipettor cones.
- Do not draw volumes of liquid exceeding the pipettor's range.
- Be sure that tips are fitted properly.
- During operation, the pipettor must be held vertically, tips should be immersed in liquid to a depth of 2 to 4 mm.

VWR Signature™ Ergonomic High-Performance Pipettor

ENGLISH

The pipetting button should be depressed and released slowly and evenly.

- When working with viscous samples, new tips must be rewetted prior to pipetting, by drawing and expelling the liquid to be measured.
- When pipetting liquids which tend to wet walls of tips, such as serums, proteins, organic solvents, measuring must be performed much more slowly than with other liquids.
- Tips must be replaced with new ones, when changing from one liquid to another, or if drops of liquid remain inside the tips.
- After work, the pipettor should be stored vertically in holder with tips removed.

5 - RECALIBRATION

VWR Signature™ Pipettors are calibrated by a gravimetric method, using **VWR** tips and distilled water, at the temperature $20 \pm 1^\circ\text{C}$, based on EN ISO 8655.

If during pipettor operation you find that the accuracy error (the difference between the real aspirated volume and the preset volume) exceeds the permissible value given in the table in section 1, the pipettor recalibration procedure should be carried out.

Before starting the recalibration procedure it is necessary to check whether the following requirements have been fulfilled during error determination:

- the ambient temperature, and the temperature of the pipettor, tips and water was identical
- the density of the liquid used is close to that of distilled water
- the balance with appropriate sensitivity has been used

Volume checked [μl]	Balance sensitivity [mg]
0.1 - 10	≤ 0.001
10 - 100	≤ 0.01
> 100	≤ 0.1

- $\text{mg}/\mu\text{l}$ conversion factor has been taken into account
- the requirements given in sections 3 and 4 have been fulfilled

If the above conditions are satisfied and the accuracy error for selected volume given in section 1 exceeds the permissible value, the pipettor recalibration procedure should be carried out.

The recalibration can be performed within one full turn of the key to the right or to the left only.

Recalibration conditions:

- Ambient temperature and the temperature of the pipettor, tips and liquid should be within the range $20\text{--}25^\circ\text{C}$ and stabilized during weighing within $\pm 0.5^\circ\text{C}$
- Measurements should be conducted using distilled water
- Balance sensitivity should be suitable for the volume to be controlled

Recalibration procedure:

- Set the dose volume depending on the pipettor volume according to the following table:

Model	Range of the pipette volumes [μl]	Preset volume [μl]	Permissible volumes [μl]	Volume change ΔV for full turn of the calibration key [μl] (24 increments)
VWR ME10 VWR MT10	1 - 10	1	0.92 - 1.08	0.33
VWR ME50 VWR MT50	5 - 50	5	4.8 - 5.2	1.67
VWR ME200 VWR MT200	20 - 200	20	19.4 - 20.6	6.30
VWR ME300 VWR MT300	50 - 300	50	49.2 - 50.8	10.00

- Perform three aspiration series (each series should include the aspirations from all channels), weigh each time and calculate the average value of the aspirations.
- Calculate average aspirated volume in μl multiplying the average aspiration amount [mg] by the distilled water density coefficient [$\mu\text{l}/\text{mg}$], which depends on temperature and pressure according to the following table:

Temperature [$^\circ\text{C}$]	Pressure [kPa]		
	95.0	101.3	105.0
20	1.0028	1.0029	1.0029
21	1.0030	1.0031	1.0031
22	1.0032	1.0033	1.0033
23	1.0034	1.0035	1.0036
24	1.0037	1.0038	1.0038
25	1.0039	1.0040	1.0040

VWR Signature™ Ergonomic High-Performance Pipettor

ENGLISH

If the average aspirated volume exceeds the permissible value, the following should be done:

- Remove the pipetting pushbutton, (Fig. 4A)
Warning: The pipetting pushbutton consists of 2 parts: the knob (Fig. 1A2) and the pushbutton (Fig. 1A1). After removal of pushbutton, both parts are separated.
- Holding the volume adjustment knob to protect it against rotation, insert the calibration key into the cuts of the calibration screw, (Fig. 4B)
- Turn the key clockwise to reduce the aspirated volume, or counter-clockwise to increase the volume, (Fig. 4C). One full turn of the calibration key changes the pipettor aspiration volume by the amount given in the table
- Take out the key and fix the pipetting push-button, (Fig. 4D). The pipetting pushbutton should be fixed in opposite order.

Determine the average aspirated volume. The average volume should be within the permissible range given in the table. If the volume exceeds the values stated, the recalibration procedure should be repeated.

6 - CLEANING AND STERILIZATION

Cleaning:

External surfaces of the pipetting pushbutton, the ejector pushbutton, the handgrip, the shafts and the adjustment knob may be cleaned using a cloth dampened in isopropyl alcohol.

Sterilization:

The pipettors can be sterilized in the autoclave at 121°C for 20 minutes maximum. Exceeding the recommended autoclaving temperature or time will cause damage to the pipettor and will void the warranty. After sterilization, the pipettor should be dried and cooled down to room temperature.

It is recommended:

- to sterilize the pipettors in autoclave with an initial vacuum and drying cycle.
- prior to sterilization to unscrew the nut connecting the handle and module slightly. After autoclaving these parts should be screwed tight again.
- to set the locking ring in lower (unlocked) position prior to sterilization.
- to test pipettor calibration every 10 sterilization cycles.

7

7 - PIPETTOR KIT AND ACCESSORIES

Pipettor kit

The following items are included in each package:

- Pipettor
- Instruction manual
- Calibration key
- Identification labels

Accessories

Modell	Tip Choice	Cat. No. US	Cat. No. EU
VWR ME10 VWR MT10	10 µl	53509-130	613-0334
VWR ME50, VWR MT50, VWR ME200 VWR MT200	200 µl	53508-783	613-0241
VWR ME300 VWR MT300	300 µl	53509-126	613-0266

8 - SPARE PARTS

The spare parts are indicated in Fig.1, 4, 5:

A: Pushbutton A1: Pushbutton A2: Knob

B: Adjustment knob

D: Ejector

E: Piston assembly

F: Ejector pushbutton

G: Calibration key

H: Shaft

P: Locking ring

Spare parts can be ordered from a **VWR International** representative (type of the pipettor and name of the parts for this pipettor should be specified). Damaged pipettors should be returned for repair to service company's facility with accompanying information on the kind of liquids that have been measured.

Warning: Replacing the piston will require conducting of calibration procedure according to section 5.

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VWR International, Inc.