

Labpette R™ Repeating Pipette and Combi-Syringes

The Labpette R Repeating Pipette provides accurate and precise repetitive pipetting from the single filling of a disposable tip. Use of the Labpette R saves time and reduces the fatigue associated with repetitive pipetting using a standard pipette.

Twenty-nine different volumes between 1 µL and 5,000 µL may be dispensed with this pipette. As many as 48 pipetting steps may be carried out without refilling the reservoir tip. Samples can be dispensed at intervals as fast as one second without compromising accuracy or reproducibility.

The slender body of the Labpette R fits comfortably in the natural contour of your hand. The plunger is located on the top of the unit, like that of a standard pipette, for easy access and reduction of hand strain.

The Labpette R is simple to use – select the appropriate tip and lock in place. Then choose the desired volume setting with the click stop adjustment dial and fill the tip with the required volume, using the sliding fill lever. To dispense, depress the plunger. For convenience, a chart showing the volume dispensed for each tip size and dial position is located on the body of the pipette.

The Labpette R performs best when used with Labnet's BioFree® Combi-Syringe™ dispenser tips. Available sterile or nonsterile, Combi-Syringes are supplied in sturdy, stackable dispenser boxes of 100.



CAT NO.		DESCRIPTION			
P3000		Labpette R Repeating Pipette, 1 µl to 5,000 µl			
BIOFREE COMBI-SYRINGE DISPENSER TIPS					
CAT NO.	CAT NO.	TIP			QUANTITY
Nonsterile	Sterile	Capacity	Accuracy	Precision	Per Box
P3505 ¹	P3505-S	0.05 mL	±2.5%	<5.0% to <2.5%	100
P3512 ²	P3512-S	0.5 mL	±0.8%	<0.7% to <0.4%	100
P3515 ¹	P3515-S	1.25 mL	±0.8%	<0.5% to <0.2%	100
P3518 ¹	P3518-S	2.5 mL	±0.7%	<0.4% to <0.2%	100
P3521	P3521-S	5.0 mL	±0.4%	<0.3% to <0.2%	100
P3524	P3524-S	12.5 mL	±0.3%	<0.3% to <0.1%	100
P3525*	P3525-S	25.0 mL	±0.3%	<0.2% to <0.1%	100
P3527*	P3527-S	50.0 mL	±0.2%	<0.2% to <0.1%	100
P3535	Combination pack, 20 each of 0.5, 1.25, 2.5, 5.0 and 12.5 mL				
P3527-A	Adapter for 25 and 50 mL Combi-Syringes				

* Four adapters included

¹ result using a 10µl tip attached

² result using a 200 µl tip attached