[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Dispensette® S Trace Analysis

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Table of contents

	13.1	Sending for repair	78
14		Calibration service	80
15		Warranty	81
16		Disposal	82
1		Introduction	44
1.1		Scope of delivery	44
1.2		Use operating manual	44
2		Safety regulations	46
2.1		General safety regulations	46
2.2		Function	47
2.3		Limitations of Use	47
2.4		Operating Limitations	47
2.5		Operating exclusions	48
2.6		Storage conditions	48
2.7		Recommended application range ..	48
3		Functional and operational compo- nents	50
4		Assembly	51
4.1		First Steps	51
4.2		Priming	53
5		Operation	56
5.1		Dispensing	56
5.2		Accessories	57
5.3		Bottle stand	60
6		Error limits	61
7		Checking the Volume (Calibration)	62
8		Adjustment	63
8.1		Making adjustments	63
9		Cleaning	64
9.1		Standard cleaning	64
9.2		Cleaning for trace analysis	65
9.3		Replacing the dosing cannula/ valves	66
9.4		Replacing the dispensing cartridge	68
10		Troubleshooting	70
11		Product markings	72
12		Ordering Information	73
13		Repairs	78

1 Introduction

1.1 Scope of delivery

Dispensette® S Trace Analysis bottle-top dispenser, for GL 45 threaded bottles, discharge tube or discharge tube with recirculation valve, telescopic aspiration tube, recirculation tube (optional for devices with back-dosing valve), assembly wrench, various bottle adapters, a quality certificate and these instructions for use.

Nominal volume ml	Adapter for bottle thread	Filling tube length mm
10	GL 28/ S 28 (ETFE), GL 32 (ETFE), GL 38 (ETFE), S 40 (PTFE)	125- 240

1.2 Use operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the instructions for use at www.brand.de/om.

1.2.1 Signal words and their meaning

Signal word	Meaning
 WARNING or  WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
 CAUTION or  CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2.2 Symbols used in the instructions

The sign warns of a possible hazard.

The sign warns of hazardous media and the dangers that can arise from liquid splashes.

The symbol indicates that the instruction manual and safety information must be observed.

The sign suggests that eye protection be worn.

The sign suggests that you to wear protective clothing.

The sign suggests that protective gloves be worn.

The sign suggests that a face shield be — — worn during work.

1.2.3 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The laboratory device Disensette® S Trace Analysis Can be used in combination with hazardous materials, operations and equipment. However, the instructions for use cannot identify all the safety problems that may occur. It is the responsibility of the user to ensure compliance with safety and health regulations and to determine the appropriate restrictions before use.

- Each user must have read and observe the instructions for use accompanying the laboratory device before using the device. The laboratory equipment must only be used by trained and qualified personnel.
- Follow general hazard warnings and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous materials), the general safety rules in the laboratory must be followed and regulations for handling the samples/media must be observed. The information provided by the media manufacturers (e.g. safety data sheets) must be observed.
- Use the laboratory equipment only for pipetting or dispensing media within the defined limits and limitations of use. Observe usage exclusions.
- When working with flammable media, take precautions to avoid electrostatic charges, e.g. do not dose into plastic containers and do not rub equipment with a dry cloth. Do not use the laboratory equipment in explosive atmospheres. If in doubt, it is essential to contact the manufacturer or dealer.
- Always check that the laboratory equipment is in good condition before use. If there are any signs of malfunctions in the laboratory equipment (e.g. Stiff piston, leaks or power supply), immediately stop working with the appliance and refer to the troubleshooting section in the user manual. If necessary, contact the manufacturer.
- Always work in such a way that neither users nor other persons are at risk. Avoid splashes. Use only suitable vessels. Never use unnecessary force or force while operating, cleaning, or maintaining the laboratory equipment.
- If the laboratory equipment is powered by a power supply, batteries or rechargeable batteries, the correct condition of the components and the connection to the equipment must be checked regularly. Do not operate the laboratory equipment or its accessories in an unprotected, damp or wet environment.
- Do not make any technical changes. Use only original manufacturer's spare parts, and do not use power supplies or batteries of identical sizes and specifications from other manufacturers. Do not disassemble the laboratory equipment and its accessories (e.g., power supplies, cables, stands, batteries, etc.) any further than the instructions for use.
- Do not autoclave the laboratory equipment unless permitted by the instructions for use.
- Autoclaving the instrument Disensette® S Trace Analysis is not permitted!

- Never press down the piston when the discharge tube is sealed with the screw cap.
- Never remove the discharge tube while the dispensing cylinder is filled.
- Reagents can accumulate in the screw cap of the discharge tube. Therefore, clean the screw cap regularly.
- To prevent tipping, use a bottle stand – particularly with small bottles and when using the flexible discharge tube.
- Never carry the instrument mounted on the bottle by the cylinder sleeve or the valve block. Breakage and detachment of the cylinder can lead to injuries from chemicals, among other things; see First Steps, p. 51, Fig. 4.
- Never use force. Always pull the piston up and press down gently when dispensing.
- Always tighten the retaining ring between the valve block and dispensing cylinder by hand. Do not use any tools.

2.2 Function

With the Dispensette® S Trace Analysis bottle-top dispenser, liquids can be dispensed directly from the supply bottle. The instruments are marked DE-M and optionally equipped with recirculation valve.

2.2.1 Handling

If handled properly, the dispensed liquid comes into contact with the following chemical-resistant materials only:

Various fluoroplastics (e.g. ETFE, FEP, PFA, PCTFE, PTFE), Al₂O₃-sapphire, platinum-iridium or tantalum depending on the version (see label on the dispensing tube and the recirculation valve).

2.3 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- operating temperature from +15°C to +40°C (from 59°F to 104°F) of instrument and reagent
- vapor pressure up to 600 mbar. Absorb slowly above 300 mbar in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm²/s (dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Density up to 3.8 g/cm³

2.4 Operating Limitations

- Liquids that form deposits (e.g., crystallizing solutions or concentrated alkaline solutions) may make the piston difficult to move or may cause jamming. If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately. See also Cleaning, p. 64.
- When dispensing inflammable media, make sure to avoid to buildup of static charge (e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth).

- If hazardous media (e.g., hydrofluoric acid or bromine) are dispensed, replace the dispensing unit after approx. 3000 complete strokes. Replacement may also be necessary earlier depending on the medium dispensed and the frequency of use. See also Replacing the dispensing cartridge, p. 68.
- The instrument is designed for general laboratory applications and complies with the requirements of the relevant standards, e.g. DIN EN ISO 8655. The use of the instrument for special applications (e.g. in trace analysis, in the food sector, etc.) must be carefully checked by the user. Approvals for specific applications (e.g., for the production and administration of food, pharmaceuticals, or cosmetics) have not been given.

2.5 Operating exclusions

* Dissolved sodium azide is permitted up to a concentration of 0.1%.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps.

2.5.1 Dispensette® S Trace Analysis

Never use Dispensette® S Trace Analysis for:

- Liquids that attack Al_2O_3 -sapphire or fluoroplastics such as ETFE, FEP, PFA, PCTFE, and PTFE (e.g., dissolved sodium azide*)
- Liquids that are decomposed catalytically by platinum-iridium (e.g., H_2O_2) or attack tantalum. Take notice of the valve spring material of the instrument version in use
- Organic solvents
- Trifluoroacetic acid
- Explosive liquids (e.g., carbon disulfide)
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument

* Dissolved sodium azide is permitted up to a concentration of 0.1%.

2.6 Storage conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

2.7 Recommended application range

Select valves with valve springs made of platinum-iridium or tantalum according to the intended use. The instrument can be used for the following dispensing media:

Dispensing medium	Valve springs: Pt-Ir	Valve springs: Ta
Acetic acid	✓	✓

Dispensing medium	Valve springs: Pt-Ir	Valve springs: Ta
Ammonia solution	✓	✓
Bromine	✓	✓
Hydrochloric acid	✓	✓
Hydrofluoric acid*	✓	—
Hydrogen peroxide	—	✓
Nitric acid	✓	✓
Perchloric acid	✓	✓
Phosphoric acid	✓	✓
Sodium hydroxide, 30%	✓	—
Sulphuric acid	✓	✓
Water	✓	✓

Legend:

✓ = The device is suitable for the medium

— = The device is not suitable for the medium.


This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

* Hydrofluoric acid slightly attacks the sapphire components (99.99% Al_2O_3) used in the valves and can release aluminum ions. Depending on the analysis method used (inorganic trace analysis), this can lead to slightly increased aluminum blank values. To reduce the aluminum values, we recommend discarding three to five doses of 2 ml each before the analysis.

Compounds that contain fluoride (e.g., NaF) attack tantalum.

Last updated: 0419/2

3 Functional and operational components

4 Assembly

⚠ WARNING

Risk of injury from hazardous media

➤ Wear your personal protective equipment, in particular protective clothing suitable for the medium to be used, eye protection and protective gloves!

➤ Before use, find out how you can use the medium safely and whether the laboratory device is suitable for it.

➤ Observe the chemical-specific safety data sheets and follow the operating instructions for the medium at the workplace!

4.1 First Steps

NOTICE Select correct discharge valve and discharge tube

Discharge valve and discharge tube are marked with the spring material. The label 'Pt-Ir' or 'Ta' must be the same on each instrument. The spring material is determined by the application range; see Recommended application range, p. 48.


When labels do not match, the desired application is not possible as components can become corroded or decomposed.

1. Check retaining ring

- a. Check that the retaining ring is screwed tight.

2. Installing the filling tube/recirculation tube

- a. Adjust the length of the telescopic filling tube to the bottle height and attach it.
- NOTICE** Set the filling tube to approx. 2 cm above the bottom of the bottle so that no crystallized substances or solids are sucked into the instrument.
- b. Center the filling tube (side with smaller diameter) and attach it carefully to avoid damaging the olive-shaped nozzle.
 - ⇒ If a dispensing cannula with recirculation valve is used, the recirculation tube must also be fitted.
 - c. Insert the recirculation tube with the opening pointing outward.

3. Mounting and aligning the instrument on a bottle

- a. Screw the instrument (GL 45 threads) onto the reagent bottle and then align the dispensing cannula with the bottle label. This is done by rotating the valve block with the dispensing cannula.

NOTICE Do not tip

To prevent tipping, use a bottle stand – particularly with small bottles and when using the flexible discharge tube.

- NOTICE Clean instrument before use in trace analysis:**
See Cleaning for trace analysis, p. 65.

4. Transporting the instrument

⚠ WARNING! Hazardous media:

If hazardous media are used, handle the instrument and bottle only with protective gloves.

⚠ WARNING! Risk of injury:

Carry the instrument mounted on the reagent bottle only as shown in the illustration. Otherwise, the bottle neck, the retaining ring, or the bottle adapter may break off.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. When mounted to a reagent bottle, always carry the instrument as shown in the figure.

NOTICE Use the correct bottle adapters

The scope of delivery includes bottle adapters made of polypropylene (PP). These may be used only for media that do not attack PP. Alternatively, ETFE/PTFE bottle adapters can be used (Accessories/spare parts). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters.

NOTICE Template for adapter

An adapter template is available for download at www.brand.de in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop <https://shop.brand.de> in the “Selection Charts” section.


4.2 Priming

NOTICE Rinse instrument thoroughly:


Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Dose slowly in order to prevent spurts. Perform cleaning for trace analysis depending on requirements (Cleaning for trace analysis, p. 65).

Instruments with recirculation valve

- a. Open the screw cap of the dispensing tube. To avoid splashes, hold the discharge tube orifice on the inner wall of a suitable receiving vessel.

- b. Turn the valve to 'Re-dose'!

- c. For priming, pull the piston up approx. 30 mm and push it down until the lower stop. Repeat this procedure at least 5 times.

- d. Turn the valve to 'Dispense'!

- e. To avoid splashing, hold the opening of the dispensing cannula against the inside of a suitable collection container and dispense until the air has been removed from the dispensing cannula and no bubbles remain. Wipe away any remaining drops from the cannula.

Instruments without recirculation valve

- a. Open the screw cap of the dispensing cannula (see 'Instruments with recirculation valve', Fig. a). To avoid splashes, hold dispensing cannula orifice on the inner wall of a suitable receiving vessel.
- b. For priming pull the piston up approx. 30 mm and push it down until the lower stop. Repeat this process about 5 times until the air has been removed from the dispensing cannula and no bubbles remain.

5 Operation

⚠ WARNING


To be observed for every application, especially with hazardous media

- Wear your personal protective equipment, in particular protective clothing, eye protection and protective gloves!
- Before use, find out how to use the medium safely and whether the laboratory instrument is suitable for this purpose. Observe the safety data sheets and follow the operating instructions for the medium!
- Never press down the piston when the dispensing cannula is sealed with the screw cap. Dispense slowly to avoid splashing.
- Open the screw cap slowly to avoid splashing. Media residue can accumulate in the screw cap.

5.1 Dispensing

1. Select volume


- a. Loosen the volume selector thumb screw $\frac{3}{4}$ turn (1), set the pointer to the desired volume (2) and then retighten the volume thumb screw (3).

2. Dispensing

- a. Unscrew the screw cap of the dispensing cannula.
- b. For instruments with a recirculation valve, turn the valve to dispensing.
- c. Hold the dispensing cannula orifice on the inner wall of a suitable receiving vessel.

- d. Gently pull the piston up to the stop and then push it back down to the lower stop evenly and without applying too much force.


- e. Wipe off the dispensing cannula against the inner wall of the receiving vessel.
- f. Seal the dispensing tube with the screw cap.

NOTICE Park position:

After use, always press the plunger down to the bottom stop (**park position**). If the piston is not pressed down to the bottom stop, media may escape unintentionally.

NOTICE Label the instrument when it is filled.

The filled state of the instrument during cleaning must be specially marked!

5.2 Accessories

5.2.1 Flexible discharge tube with recirculation valve

⚠ WARNING

Risk of media contact if the hose is damaged

If the hose is damaged, media contact may occur. Check the hose for damage (e.g. kinks and the like) before each use and replace it if necessary.

- Use a face shield and personal protective equipment.
- Secure the bottle with a bottle stand.
- In order to avoid splashing reagent, hold the dispensing tube firmly, and insert it into the holder provided after use.
- Rinse the tubing to clean it.
- Do not disassemble.


With the exception of HF, the flexible discharge tube can be used for serial dispensing (Accessories/ spare parts).

The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, make sure that the hose is properly looped and not twisted. The operating exclusions for the respective instrument in use apply.

⚠ WARNING! Not suitable for hydrofluoric acid (HF)


The flexible discharge tubing must not be used for dispensing HF (hydrofluoric acid)!

NOTICE Pay attention to correct valve marking see see Replacing the discharge valve, p. 67.


Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, p. 64).
- a. For instruments with a recirculation valve, set the valve to 'Recirculation' and pull the valve toggle upwards.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.

- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it. The Instrument must not be mounted on a bottle for this purpose. Fit the collecting tube.
- d. Press down the plug of the recirculation valve.

- e. Slide the flexible discharge tube housing into the valve block up to the stop.

- f. Push the housing all the way down.
- g. Attach the valve handle that fits the discharge valve and press it in firmly. Note the color coding and marking.

NOTICE Use bottle stand, see accessories/spare parts.

5.2.2 Drying tube

For moisture- or CO₂-sensitive media, it may be necessary to use a dry pipe filled with suitable absorbent (not included in the scope of delivery).

See also Accessories/spare parts.

- a. Use a coin to unscrew the air vent cap.

- b. Screw in the filled drying tube.
- c. Place the PTFE sealing ring on the bottle thread or place the screwed-on bottle adapter and screw the Instrument onto the bottle.

NOTICE If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

5.2.3 Sealing ring for valve block

For highly volatile or outgassing media, we recommend sealing the connection from the valve block to the bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts).

5.3 Bottle stand

For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over (Accessories/spare parts).

Assembly

- a. Position the mounting plate at the appropriate height.
- b. As shown, insert the instrument firmly into the holder until it audibly clicks into place.
- c. Then lock the holder into place with the screw.

6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. The test was carried out in accordance with DIN EN ISO 8655-6 with a completely filled instrument and uniform and jerk-free dispensing.

Error limits

Nominal volume ml	A* ≤ ± %	µl	CV* ≤ %	µl
10	0.5	50	0.1	10

*A = Accuracy, CV = Coefficient of Variation

Partial volumes

The percentage values for A and CV are relative to the nominal volume (V_N) and must be converted for partial volumes (V_P).

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

e.g.	Volume	A* ≤ ± %	µl	CV* ≤ %	µl
V_N	10.0	0.5	50	0.1	10
$V_T = 50\% N$	5.0	1.0	50	0.2	10
$V_T = 10\% N$	1.0	5.0	50	1.0	10

*A = Accuracy, CV = Coefficient of Variation

NOTICE

The error limits in DIN EN ISO 8655- 5 are significantly lower. The maximum error for a single measurement can be approximated from the sum of error limits $EL = A + 2 \times CV$ (for 10 ml: $50 \mu\text{l} + 2 \times 10 \mu\text{l} = 70 \mu\text{l}$).

7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3–12 months. This time frame should be adjusted to correspond with individual requirements.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits, p. 61) is performed as follows:

1. Prepare instrument

Clean the instrument (Cleaning, p. 64), fill it with test liquid (distilled water) and then prime it carefully.

2. Check volume

- Carry out 10 doses with test liquid in 3 volume ranges (100 %, 50 %, 10 %).
- For discharge, depress the piston slowly and steadily without force until the lower stop
- Wipe off the tip of the dispensing tube.
- Weigh the dispensed quantity on an analytical balance. Refer to the operating manual of the balance manufacturer.
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = Weighing results

N = number of weighings

V_0 = nominal volume

Z = correction factor (e.g. 1.0029 $\mu\text{l}/\text{mg}$ at 20 °C, 1013 hPa)

Mean value:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$R\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$VK\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

8 Adjustment

After a long period of usage, an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (Checking the Volume (Calibration)).
- Calculate mean volume (actual value) (Checking the Volume (Calibration)).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm the appropriate adjustment.

8.1 Making adjustments

- Insert the pin of the installation tool into the calibration cover (pos. 1) and break it off by turning it. Dispose of the calibration cover.

- Insert the pin of the installation tool into the calibration screw and turn to the left to increase the dispensing volume or turn to the right to decrease the dispensing volume (e.g. current value 9.97 ml approx. 1/2 turn to the left).

Nominal volume	Analog/Fixed max. +/- [μ l]	One rotation corresponds to [μ l]
10 ml	60	~ 80

- The calibration is complete.
- ⇒ The change in calibration is indicated by a red disk (circle in illustration).

9 Cleaning

⚠ WARNING

There is a risk of exposure to liquids during cleaning.

The components cylinder, valves, suction pipe and dosing cannula are filled with liquid!

➤ Wear protective clothing, eye protection and protective gloves to protect yourself from hazards caused by the liquid during cleaning.

➤ Never point the openings of the suction tube, dosing cannula and valves towards the body in order to avoid hazards from escaping liquid.

➤ Never remove the dosing cannula when the dosing cylinder is full, as liquid will then escape.

The instrument must be cleaned in the following situations to assure correct operation:

- Before first use
- Immediately if the piston becomes stiff
- Prior to a reagent change
- Prior to long term storage
- Prior to dismantling the instrument
- Before valve replacement
- Regularly when using liquids which form deposits (e.g., crystallizing liquids)
- Regularly if liquid has accumulated in the screw cap

The instrument must **not** be autoclaved!

9.1 Standard cleaning

1. Empty the instrument completely

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing.

NOTICE Recirculation

If the instrument is equipped with a recirculation valve, it must be emptied in dispensing and recirculation positions.

2. Rinse the instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and fill and empty completely several times to rinse.

NOTICE Recirculation

If a recirculation valve is present, also fill the instrument completely several times in the recirculation position and empty via recirculation.

3. Rinsing in recirculating position (optional)

If the instrument is equipped with a recirculation valve, it must also be rinsed in recirculating position after rinsing the instrument.

- a. Set the recirculation valve to "recirculate" and fill and discharge the instrument several times.

9.2 Cleaning for trace analysis

Before use in trace analysis, the instrument must first be thoroughly cleaned. In doing so, use reagents having the purity level "per analysis" or better. If contamination of the bottle contents should be avoided, do not use the instrument with a recirculation valve. If the instrument is used with a recirculation valve, the cleaning must be done using both the dispensing and recirculating function (see Figures 1 and 2).

Good results have been achieved in practice with the following recommended cleaning methods. Modify accordingly, as necessary.

- a. Screw the device onto a bottle filled with **acetone**, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap. Dispense twice after approx. 24 hours reaction time, then the device must be completely drained and rinsed 5 times with pure water.
- b. Screw the device onto a bottle filled with approx. 20 % **hydrochloric acid**, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap.
- c. After approx. 24 hours reaction time, dispense twice and fill again to the maximum.
- d. Repeat step 3 two more times. After another 24 hours of reaction time, completely drain the device and rinse it 5 times with pure water.
- e. Repeat steps 2 to 4 with a solution of approx. 30% nitric acid.
- f. Screw the device onto a bottle filled with the desired dispensing medium, deaerate and fill up to the maximum. Leave the piston at the upper stop and close the discharge tube with the screw cap.
- g. After approx. 24 hours reaction time, dispense twice and fill again to the maximum.

1 Dosing function

2 Recirculation function

- h. Repeat step 7 two more times. After an additional 24 hours reaction time, dispense twice and leave the piston at the lower stop.

NOTICE If cleaning is insufficient, repeat the cleaning procedure.

9.3 Replacing the dosing cannula/valves

⚠ WARNING

The valves, telescoping filling tube and discharge tube contain reagent!

Contact with potentially hazardous media is possible.


- Clean the instrument before the discharge tube, valves or dosing unit are exchanged.
- Do not disassemble the dosing unit.
- Follow the safety instructions (Safety regulations, p. 46).

NOTICE A functional test must always be carried out after replacing components.

9.3.1 Replacing the dosing cannula

- a. For instruments with a recirculation valve, set the valve to 'Recirculation'.
- b. Pull the valve toggle upwards.

- c. Slide the dispensing cannula housing all the way up, then pull it forward with gentle up and down motions.

- d. Slightly push the dosing cannula onto the valve.

Then pull the dosing cannula slightly upwards.
The coupling piece slides downwards

Push the dosing cannula onto the valve as far as it will go.

Press the dosing cannula down to lock it.

- e. Push the housing all the the discharge tube way down.
- f. For units with a redosing valve, place the valve knob in the 'redosing' position and push it downwards.
- g. Carry out a functional test after replacement.

9.3.2 Replacing the valves

9.3.2.1 Replacing the discharge valve

- a. After disassembling the dispensing cannula (see Replacing the dosing cannula, p. 66), use the assembly tool to unscrew the discharge valve.

- b. Firmly screw in the new discharge valve first by hand and then tighten it with the assembly tool. The thread should not be visible.

NOTICE Always install valves intended for the respective instrument type

The filling valve of the Disensette® S Trace Analysis is the same for all instrument versions, however, the discharge valve is different. Make sure to use only the filling valve marked with 'S' (sapphire). For easy identification, the discharge valves for the Disensette® S Trace Analysis are marked with 'Pt-Ir' or 'Ta'.

9.3.2.2 Replacing the filling valve

- a. Pull out the recirculation tube and the telescopic filling tube.

- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

9.3.2.3 Loosen stuck filling valve

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 µl plastic pipette tip (see the figure at the side).

9.4 Replacing the dispensing cartridge

The dispensing unit is a wear part and must be replaced according to the specific conditions of use. The replacement interval depends on the medium used and number of dispensing cycles, whereby replacement is to be expected at the latest after 10,000 full strokes (medium: distilled water). Please note that media residues can accumulate inside the dispensing unit; these can be released by tilting the instruments or dispensing unit. If drops of media are observed on the outside of the dispensing cylinder or in the retaining ring, immediately check the instrument and replace the dispensing unit. If hazardous media (e.g., hydrofluoric acid or bromine) are dispensed, replace the dispensing unit after approx. 3000 complete strokes. Replacement may also be necessary earlier depending on the medium dispensed and the frequency of use.

⚠ WARNING

The valves, telescoping filling tube and discharge tube contain reagent!

Contact with potentially hazardous media is possible.

- Clean the instrument before the discharge tube, valves or dosing unit are exchanged.
- Do not disassemble the dosing unit.
- Follow the safety instructions (Safety regulations, p. 46).

NOTICE A functional test must always be carried out after replacing components.

- a. Completely unscrew the retaining ring of the dispensing unit by hand and remove the dispensing cartridge.
- b. Insert a new dispensing unit into the valve block and tightly screw on the retaining ring by hand.
- c. Perform a functional test for leaks.

NOTICE

Factory adjustment

The dispensing unit is pre-adjusted in the factory. Therefore, calibration is not necessary after replacement.

- Before using the instrument, rinse thoroughly and discard the first dispensed volume.
- Avoid splattering.
- Depending on the requirement, perform cleaning for trace analysis (Cleaning for trace analysis, p. 65).
- Disassembly of the dispensing unit is not permitted.


10 Troubleshooting

Problem	Possible cause	Corrective action
Liquid is above the piston	Piston leaky	Perform standard cleaning (Standard cleaning, p. 64), replace dispensing cartridge (Replacing the dispensing cartridge, p. 68).
Piston sluggish	Crystal deposits, contamination	Stop dispensing immediately. Perform standard cleaning (Standard cleaning, p. 64).
Filling not possible	Volume setting screw at the lower stop	Set the desired volume (Dispensing, p. 56).
	Filling valve stuck	Clean the filling valve. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Loosen stuck filling valve, p. 68). If necessary, replace the filling valve.
Dispensing not possible	Discharge valve stuck	Unscrew the discharge valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it. If necessary, replace the filling valve (Replacing the discharge valve, p. 67).
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently	Discharge valve is not screwed in deeply enough	Tighten the discharge valve with the assembly tool until it meets the stop so that the threads are no longer visible.
Air bubbles in the instrument	Reagent with high vapor pressure has been drawn in too quickly	Slowly draw in reagent.
	Screw connections loose	Tighten the valves firmly with the assembly tool
	The instrument has not been primed	Prime the instrument (Priming, p. 53).
	Filling tube is loose or damaged	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged	Perform cleaning (Cleaning, p. 64). Tighten the valves firmly with the assembly tool.
	Recirculation tube not connected	Install the recirculation tube (First Steps, p. 51).
Dispensed volume too low	Filling tube is loose or damaged	Perform cleaning (Cleaning, p. 64). Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.

Problem	Possible cause	Corrective action
	Filling valve contaminated, loose or damaged	Perform cleaning (Cleaning, p. 64). Tighten the valves using the assembly tool. If necessary, replace filling valves.
Leaking liquid at the retaining ring	Dispensing cartridge loose or piston seal damaged	Tighten the retaining ring. If necessary, replace the dispensing cartridge
Leaking liquid between instrument and bottle	Recirculation tube not connected	Install the recirculation tube (First Steps, p. 51).
	Volatile reagent dispensed without sealing ring	Install the seal (Accessories, p. 57)

11 Product markings

Symbol or number	Meaning
	General warning symbol
	Follow the operating instructions
	Use eye protection
	Use hand protection
	Wear protective clothing
XXZXXXXX	Serial number
 25	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
www.brand.de/ip	Patent information

12 Ordering Information

Dispensette® S Trace Analysis, Analog

Volume ml	Valve spring	Without recirculation valve Art. No.	With recirculation valve Art. No.
10	Platinum-iridium	4640040	4640041
10	Tantalum	4640240	4640241

See Scope of delivery, p. 44.

Bottle adapter

External thread	for bottle thread/ground joint	Material	Art. No.
GL 32	GL 25	ETFE	704375
GL 32	GL 28/ S 28	ETFE	704378
GL 32	GL 45	ETFE	704395
GL 45	GL 32	ETFE	704398
GL 45	GL 38	ETFE	704399
GL 45	S* 40	PTFE	704391

* buttress thread

NOTICE Template for adapter:

An adapter template is available for download at www.brand.de in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop store.brand.de in the "Selection Charts" section.

Discharge tubes

With and without recirculation valve.

Nominal volume 10 ml.

Screw cap ETFE.

Discharge tube marked with 'Pt-Ir' or 'Ta'.

Packaging unit 1 piece.

Valve spring	Length mm	Without recirculation valve Art. No.	With recirculation valve Art. No.
Platinum-iridium	105	708022	708122
Tantalum	105	708024	708124

Flexible discharge tube with recirculation valve

PTFE, spiraled, approx. 800 mm long with safety handle.

Packaged unit Pack of 1.

Not suitable for hydrofluoric acid (HF).

Nominal volume ml	Dispensing hose Outer diameter	Dispensing hose Inner diameter	Art. No.
10	3	2	708132

Dispensette® S Trace Analysis discharge valve

PFA/sapphire

Valve marked with 'Pt-Ir' or 'Ta'

Packaged unit 1 pc

for nominal volume ml	Valve spring	Art. No.
10	Platinum-iridium	6732
10	Tantalum	6733

Dispensette® S Trace Analysis filling valve

PFA/sapphire

Valve marked with 'S'

Packaged unit 1 pc

for nominal volume ml	Art. No.
10	6739

**Telescopic filling tube**

FEP. Individually adjustable lengths. Packaging unit Pack of 1.

for nominal volume ml	Outer diameter mm	Length (mm)	Art. No.
10	6	70-140	708210
"	"	125-240	708212
"	"	195-350	708214
"	"	250-480	708216

**Dispensing cartridge**

Dispensing cartridge with retaining ring. Nominal volume 10 ml, adjusted, includes performance certificate.

Packaged unit 1 pc


Description	Art. No.
Dispensing unit	708035

**Venting stopper for microfilter with Luer cone**

Description	Packaging unit	Art. No.
Venting stopper for microfilter with Luer cone. PP. Venting stopper and PTFE-sealing ring.	1 pc.	704495

Recirculation tube

Description	Packaging unit	Art. No.
Recirculation tube, FEP	1 pc.	6747

Adjustment/installation tool

Description	Packaging unit	Art. No.
Adjustment/installation tool	1 pc.	6748

Sealing ring for valve block


Description	Packaging unit	Art. No.
Sealing ring for valve block. PTFE, for highly volatile media.	1 pc.	704486

Bottle stand

Description	Packaging unit	Art. No.
Bottle stand PP. Support rod, 325 mm, base plate 220 x 160 mm.	1 pc.	704275

Drying tube

Description	Packaging unit	Art. No.
Drying tube, including PTFE sealing ring (without drying agent)	1 pc.	707930

Screw cap with fastener

Packaging unit Pack of 1.

Description	Nominal volume ml	Art. No.
ETFE	10	706029

13 Repairs

13.1 Sending for repair

NOTICE Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

14 Calibration service

ISO 9001 and GLP guidelines require regular testing of your volumetric measuring devices. We recommend that you perform a volume check every 3-12 months. The cycle depends on the individual requirements of the device. In case of high frequency of use or aggressive liquids, more frequent testing should be carried out.

The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or our accredited calibration laboratory. Simply send us the device to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days. The instruments are supplied with a detailed calibration certificate or a calibration slip in accordance with DIN EN ISO/IEC 17025. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on www.brand.de (see Service and Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

15 Warranty

We shall not be liable for the consequences of improper handling, use, maintenance, operation or unauthorized repair of the device or for the consequences of normal wear and tear, in particular of wear parts such as pistons, seals, valves or glass breakage. The same applies to non-observance of the instructions for use. In particular, we assume no liability for damage caused if the device has been disassembled further than described in the user manual or if external accessories or spare parts have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

16 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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