



April, 16th 2018

Subject: Change to silicone holder for Labnet FastPette Pro and FastPette V-2 pipette controllers

Dear Valued Customer,

Thank you for continuing to choose Corning products. I would like to inform you about the change to silicone holder for Labnet FastPette Pro (Cat. No. P2002) and FastPette V-2 (Cat. No. P2000) pipette controllers as well as the silicone holder offered as spare part – Cat. No. SP29054 and SP29052.

The change will further decrease forces required for secured pipet mounting and will be implemented after stock of current products depletes.



Current



New