

Thermo Scientific Finnpiquette F3 Pipettes*

Finnpiquette F3 Variable Volume Single Channel Pipettes



Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4640000	1-10 µL	10	±0,100	0,050	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4640010	1-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,100	0,065	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4640020	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4640030	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4640080	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4640090	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4640040	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4640050	20-200 µL	200	±1,2	0,6	±1,60	0,60
		100	±1,0	0,5	±1,60	0,60
		20	±0,6	0,3	±1,60	0,60
4640060	100-1000 µL	1000	±5,0	2,0	±8,0	3,0
		500	±5,0	2,0	±8,0	3,0
		100	±1,5	0,6	±8,0	3,0
4640100	0,5-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		500	±10,0	4,0	±40,0	15,0
4640070	1-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		1000	±20,0	8,0	±60,0	30,0

Finnpiquette F3 Fixed Volume Single Channel Pipettes



Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4650000	1 µL	1	±0,040	0,040	±0,050	0,050
4650010	5 µL	5	±0,070	0,070	±0,125	0,075
4650020	10 µL	10	±0,090	0,080	±0,120	0,080
4650090	20 µL	20	±0,14	0,10	±0,20	0,10
4650030	25 µL	25	±0,15	0,125	±0,50	0,20
4650040	50 µL	50	±0,30	0,20	±0,50	0,20
4650050	100 µL	100	±0,40	0,30	±0,80	0,30
4650100	200 µL	200	±0,80	0,60	±1,60	0,60
4650060	250 µL	250	±1,0	0,8	±4,0	1,5
4650070	500 µL	500	±1,5	1,5	±4,0	1,5
4650080	1000 µL	1000	±3,0	3,0	±8,0	3,0

* Please note: Availability may vary by country.

Finnpipette F3 Multichannel Pipettes

Cat. No.	Range	Channels	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4660010	5-50 µL	8	50	±0,75	0,35	±1,00	0,40
			25	±0,63	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4660020	30-300 µL	8	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4660040	5-50 µL	12	50	±0,75	0,35	±1,00	0,40
			25	±0,63	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4660050	30-300 µL	12	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0



Thermo Scientific Finnpipette Digital Pipettes*

Finnpipette Digital Variable Volume Single Channel Pipettes

Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4500000	0.2-2 µL, micro	2	±0,050	0,040	±0,080	0,040
		1	±0,035	0,035	±0,080	0,040
		0,2	±0,024	0,020	±0,080	0,040
4500010	0.5-10 µL, micro	10	±0,100	0,050	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4500020	0.5-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4500080	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4500100	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,125	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4500110	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4500090	20-200 µL	200	±1,2	0,4	±1,60	0,60
		100	±1,0	0,4	±1,60	0,60
		20	±0,36	0,14	±1,60	0,60
4500120	100-1000 µL	1000	±6,0	2,0	±8,0	3,0
		500	±4,0	1,5	±8,0	3,0
		100	±1,0	0,6	±8,0	3,0
4500060	1-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		1000	±8,0	3,0	±40,0	15,0
4500070	2-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		2000	±20,0	6,0	±60,0	30,0



* Please note: Availability may vary by country.