

Sarpette® M

Manual pipettes in single- and multi-channel configurations



- Exceptional comfort, accuracy, and precision
- Continuously adjustable volume from either the dispenser button or thumbwheel
- Easy-grip removable and adjustable tip ejector accommodates a variety of tips and tubes
- Ergonomic, lightweight and well-balanced; comfortable for right- or left-handed use
- Low dispensing and ejecting forces
- 8 and 12-channel models feature curved ejector bars to further reduce ejecting forces
- Insulated pipette body protects micrometer from thermal expansion effect
- Fully autoclavable for easy sterilization

Ordering Information

Sarpette® Manual Pipettes	
90.1771.001	Sarpette® M1, Single-Channel, 0.1-2µl
90.1771.002	Sarpette® M1, Single-Channel, 0.5-10µl
90.1771.003	Sarpette® M1, Single-Channel, 2-20µl
90.1771.004	Sarpette® M1, Single-Channel, 10-100µl
90.1771.005	Sarpette® M1, Single-Channel, 20-200µl
90.1771.006	Sarpette® M1, Single-Channel, 100-1000µl
90.1771.007	Sarpette® M1, Single-Channel, 1-5ml
90.1771.008	Sarpette® M1, Single-Channel, 1-10ml
90.1772.001	Sarpette® M8, Eight-Channel, 1-10µl
90.1772.002	Sarpette® M8, Eight-Channel, 5-50µl
90.1772.003	Sarpette® M8, Eight-Channel, 20-200µl
90.1773.001	Sarpette® M12, Twelve-Channel, 1-10µl
90.1773.002	Sarpette® M12, Twelve-Channel, 5-50µl
90.1773.003	Sarpette® M12, Twelve-Channel, 20-200µl



Sarpette® Stands

92.1770.070	Carousel stand for 6 Sarpette® M1 (single-channel models only)
92.1770.080	Metal linear stand for 6 Sarpette® M
92.1770.090	Shelf clip for 3 Sarpette® M

Each pipettor comes with instruction manual, QC certificate, calibration tool, shelf clip, ID labels, and Sarstedt tip compatibility insert.