



## The right pipette tip for your Transferpette®

Optimal results are achieved in combination with the BRAND micropipettes. Quality tips from BRAND provide a perfect fit for precise analysis.

Pipette and filter tips from BRAND are tested for BRAND pipettes and for most pipette models from GILSON®, Thermo Fisher Scientific FINNPIPETTE®, Eppendorf®, and Sartorius®.

The 5 ml tip can be used with the pipette models from models from BRAND and Thermo Fisher Scientific FINNPIPETTE®. The 1250 µl and 10 ml tips are suitable for BRAND and Eppendorf®.

**Note:** Pipette shafts may be subject to modifications. Compatibility depends, among other things, on the manufacturer, pipette type, serial number, and manufacturing date.

We recommend using the free sample bag to check tip compatibility.



All information about the range of tip products can be found at [shop.brand.de](https://shop.brand.de)





Transferpette® S



Transferpette® electronic

Pipette tips Filter tips Volume range	Transferpette® single-channel Nominal volume											Transferpette® multi-channel Nominal volume								
	1 µl	2,5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	500 µl	1000 µl	2 ml	5 ml	10 ml	10 µl	50 µl	100 µl	200 µl	300 µl	1250 µl	
0.1 - 20 µl	✓	✓ <sub>2)</sub>	✓ <sub>2)</sub>											✓						
0.5 - 20 µl		✓	✓											✓						
1 - 50 µl		✓ <sub>2)</sub>	✓											✓						
2 - 200 µl <sup>1)</sup>				✓	✓	✓	✓	✓							✓	✓	✓	✓		
5 - 300 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓		
50 - 1000 µl <sup>1)</sup>									✓	✓										
50 - 1250 µl <sup>3)</sup>																				✓
50 - 1250 µl XL <sup>3)</sup>																				✓ <sub>2)</sub>
0.5 - 5 ml												✓	✓							
1 - 10 ml													✓							
0.1 - 1 µl	✓	✓ <sub>2)</sub>	✓ <sub>2)</sub>											✓						
0.5 - 10 µl		✓	✓											✓						
1 - 20 µl			✓											✓						
2 - 20 µl				✓	✓	✓	✓	✓							✓	✓	✓			
5 - 100 µl				✓	✓	✓	✓	✓							✓	✓	✓	✓		
5 - 200 µl								✓							✓	✓	✓	✓		
50 - 1000 µl									✓	✓										
50 - 1250 µl <sup>3)</sup>																				✓
50 - 1250 µl XL <sup>3)</sup>																				✓ <sub>2)</sub>



✓ = Tip volume less than pipette's nominal volume

<sup>1)</sup> Bulk tips are yellow or blue colored; racked tips are colorless in a yellow-colored or blue-colored tip tray

<sup>2)</sup> Calibration and possible adjustment needed

<sup>3)</sup> 1250 µl = Length 86 mm, 1250 µl XL = Length 102 mm

**BRAND GMBH + CO KG**

P.O. Box 1155 | 97861 Wertheim | Germany

T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de



**BRAND. For lab. For life.®**

BIO-CERT®, BRAND®, BRAND. For lab. For life.®, Transferpette®, as well as the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. The BRAND-GROUP figurative mark is registered trademark or trademark of Brand Group SE & Co. KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners.

Our technical literature is intended to inform and advise our customers. However, the validity of general empirical values, and of results obtained under test conditions, for specific applications depends on many factors beyond our control. Please appreciate, therefore, that no claims can be derived from our advice. The user is responsible for checking the appropriateness of the product for any particular application.

California Residents: For more information concerning California Proposition 65, please refer to [www.brand.de/calprop65](http://www.brand.de/calprop65).

Subject to technical modification without notice. Errors excepted.



Find accessories and replacement parts, user manuals, test instructions (SOP) and product videos at [shop.brand.de](http://shop.brand.de)



Further information on products and applications can be found on our YouTube channel: [mylabBRAND](https://www.youtube.com/mylabBRAND)



**Renewable energy**

In our factory we use 100% eco-friendly energy from certified hydro-electric power plants and energy generated on premise from our high efficiency cogeneration unit.



**Sustainable packaging**

We use cardboard with approx 90% recycled content for our product packaging.



**Waste reduction**

The amount of waste is reduced as follows compared to the products mentioned: TipRack compared to TipBox by more than 20%.



994792 © 2022 BRAND GMBH + CO KG | Printed in Germany | 0324

BRAND (Shanghai) Trading Co., Ltd.  
Shanghai, China

T +86 21 6422 2318  
info@brand.com.cn  
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.  
Mumbai, India

T +91 22 42957790  
customersupport@brand.co.in  
www.brand.co.in

BRANDTECH® Scientific, Inc.  
Essex, CT. United States of America

T +1 860 767 2562  
info@brandtech.com  
www.brandtech.com