



Well within Reach

The Eppendorf Pipette Holder System



Pipette Stand 2 and Charger Stand 2



Protect your precious pipette from possible damage, contamination or chemicals spilled on the bench. Our stands for single instruments have a robust design and securely hold one mechanical or electronic device. With the Charger Stand 2, your Eppendorf Xplorer®, Xplorer plus or Repeater® E3/E3x is always fully charged and ready for use.

Our Pipette Stand 2 carries one Repeater M4 multi-dispenser. This stand is also suitable for one Eppendorf Research® 3 neo, Eppendorf Research® plus, or Eppendorf Reference® 2 pipette, when used with the corresponding holder.



Details Crafted for Your Needs

Holders are available for all current Eppendorf pipettes and handheld dispensers, as well as most of their predecessors.

Pictograms clearly label pipette holders and charger shells for easy allocation in the lab.

Quickly interchangeable pipette holders and charger shells, no tool required.



Downwards compatible holders for Eppendorf Research and Reference pipettes.



Large carrying handle for easy transport of carousels from bench to bench.



Large rubber feet protect carousels and stands from liquids spilled on lab bench.



Pipette Holder 2 and Charger Shell 2

Choose from a broad range of holders designed for all handheld manual and electronic Eppendorf liquid handling instruments, as well as most predecessors. Conveniently place your pipettes or Repeater dispensers wherever you need them: Mounted to a wall,

on the shelf above your bench or inside a biological safety cabinet. All holders are compatible with the Pipette Carousel 2, and Charger Shells for installation into Charger Carousel 2 are available too.

Robust electrical contacts designed for long operational lifetime without the need for replacement.




V-shaped charging contacts ensure stable electrical connection between charger shell and electronic device.



Cord wrap functionality for unused power cord keeps your lab bench nice and tidy.



A small footprint with a diameter of 18.5 cm makes stands and carousels ideal when space is tight.



Intelligent charging electronics keep an eye on up to 6 Eppendorf Xplorer or Repeater E3/E3x dispensers (charges 4 at a time).



Magnetic quick-release power connector protects pipettes and Repeater dispensers from falling.

Pipette Carousel 2 and Charger Carousel 2

Keep all your Eppendorf pipettes and Repeater dispensers well within reach. The Pipette Carousel 2 for mechanical pipettes and the Charger Carousel 2 for electronic devices offer more flexibility and more capacity than ever before.

The Eppendorf Pipette Carousel 2 carries up to 6 mechanical single- or multi-channel Eppendorf pipettes. Different pipette holders allow for compatibility with all current and most predecessor pipettes. This includes Eppendorf Research, Research plus, Research 3 neo, Reference, Reference 2, Repeater M4 and Biomaster®.

The Charger Carousel 2 offers 50% more capacity than before and holds up to 6 Eppendorf Xplorer/Xplorer plus pipettes. When used with optional charger shells, the carousel can also carry and charge up to 6 Repeater E3/E3x or Repeater stream/Xstream multi-dispensers.

In addition to electronic devices, the Charger Carousel 2 is suitable for mechanical Eppendorf pipettes, when used with the fitting pipette holders.

More flexibility

- > Carousels carry all kinds of Eppendorf pipettes and multi-dispensers
- > With the right holders, various Eppendorf pipettes and dispensers can be organized in a single Carousel.



More capacity

- > Charger Carousel 2 for up to 6 electronic devices
- > Carousels carry both mechanical and electronic devices

More convenience

- > Interchangeable pipette holders, no tools required
- > Cord wrap functionality for unused power cord

More user safety

- > Magnetic power connector for fast and easy connection
- > Large rubber feet protect carousels from liquids spilled on bench top



> Additional information available at:
www.eppendorf.link/pipetteholder


Ordering information

Holding systems without charging feature	Order No.
Pipette Carousel 2, white , with 6 holders for Eppendorf Research® 3 pipettes, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116 000 236
Pipette Carousel 2 , with 6 holders for Eppendorf Research®, Research plus, Research plus Move It®, Eppendorf Reference®, Reference 2, Biomaster® pipettes or epMotion® Tip Tools, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116 000 015
Pipette Stand 2 , with pipette holder for one Repeater® M4 multi-dispenser, without charging function, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116 000 058
Pipette Holder 2, white , for one Eppendorf Research® 3 pipette, Pipette Carousel 2, Charger Carousel 2, Stand 2 or wall mounting, sticky tape included	3116 000 295
Pipette Holder 2 , for one Eppendorf Research®, Research plus, Eppendorf Reference®, Reference 2, Biomaster® pipette or epMotion® Tip Tool, for Pipette Carousel 2, Charger Carousel 2, Stand 2 or wall mounting, sticky tape included	3116 000 112
Pipette Holder 2 , for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, for Pipette Carousel 2, Charger Carousel 2, Stand 2 or wall mounting, sticky tape included	3116 000 120
Pipette Holder 2 , for one Repeater® E3/E3x or Repeater® stream/Xstream electronic multi-dispenser, for Pipette Carousel 2, Charger Carousel 2, Stand 2 or wall mounting, sticky tape included	3116 000 139
Pipette Holder 2 , for one Repeater® M4 multi-dispenser, for Pipette Carousel 2, Charger Carousel 2, Stand 2 or wall mounting, sticky tape included	3116 000 147
Holding systems with charging feature	
Charger Carousel 2 , with 6 charger shells for Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipettes and mains/ power adapter, additional charger shells and holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116 000 023
Charger Stand 2 , with charger shell for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, operated with mains/power adapter supplied with Eppendorf electronic pipettes and dispensers	3116 000 031
Charger Stand 2 , with charger shell for one Repeater® E3/E3x or Repeater® stream/Xstream electronic multi-dispenser, operated with mains/power adapter supplied with Eppendorf electronic pipettes and dispensers	3116 000 040
Charger Shell 2 , for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, for Charger Carousel 2 or Charger Stand 2, with charging function	3116 602 007
Charger Shell 2 , for one Repeater® E3/E3x or Repeater® stream/Xstream electronic multi-dispenser, for Charger Carousel 2 or Charger Stand 2, with charging function	3116 603 003

Learn more about the Pipette Holder System

The Eppendorf Pipette Holder System provides highly flexible carousels, stands and wall mount devices for current manual and electronic Eppendorf pipettes and multi-dispensers and most predecessors.

> www.eppendorf.link/pipetteholder



Your local distributor: www.eppendorf.com/contact

Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com