

Excel™ Electronic Pipettors

The Excel is an ergonomic, intuitive electronic pipettor for use in most laboratory applications. Its high precision stepping motor controls piston action, resulting in exceptional accuracy and precision.

Lightweight construction and ergonomic design with comfortable finger hook further reduce the risk of repetitive stress injuries. The pipettors are powered by rechargeable Li-ion battery allowing for 2,000 pipetting cycles after fully charged.

The Excel helps to optimize every lab work with functions such as multi-dispensing, serial dilution and mixing, easily navigated through the user-friendly interface. Available in single-, 8- and, 12-channel configurations, the Excel pipettors feature adjustable tip ejector along with wide tip compatibility. The Labnet Excel multi-channel pipettors are ideal for working with a 96-well format.

Each pipettor is individually tested and supplied with a certificate of quality.

- Factory calibrated to strict ISO8655 standard
- Large LCD screen and simple interface
- Ergonomically designed handle and finger rest
- Lower portion easily removed for autoclaving
- Rechargeable Li-ion battery
- 1-year warranty



Ordering Information

Cat. No.	Volume Range (µL)	Channels	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)
P3600L-10	0.5-10	1	±4.0 to ±1.0	<2.5 to <0.4	10
P3600L-20	2-20	1	±4.0 to ±1.0	<2.5 to <0.4	10
P3600L-200	10-200	1	±2.5 to ±0.6	<1.0 to <0.15	200
P3600L-1200	100-1,200	1	±2.5 to ±0.5	<0.6 to <0.15	1,000/1,200
P3608L-10	0.5-10	8	±4.0 to ±1.5	<2.5 to <0.4	10
P3608L-20	2-20	8	±4.0 to ±1.0	<2.5 to <0.3	10
P3608L-200	10-200	8	±3.0 to ±0.6	<1.0 to <0.15	200
P3608L-1200	100-1,200	8	±3.0 to ±0.5	<0.8 to <0.2	1,000/1,200
P3612L-10	0.5-10	12	±4.0 to ±1.5	<2.5 to <0.5	10
P3612L-20	2-20	12	±4.0 to ±1.0	<2.5 to <0.5	10
P3612L-200	10-200	12	±3.0 to ±0.6	<1.0 to <0.15	200
P3612L-1200	100-1,200	12	±3.0 to ±0.5	<0.8 to <0.2	1,000/1,200

Pipettors are delivered with universal power supply and set of the plugs: US, EU, UK, and AU.

Accessories

Cat. No.	Description
P3600L-BAT-1	Extra Li-ion battery, for all voltages
P3628L	Shelf clip for Excel Electronic Pipettor, holds 1 pipettor
P3630L	Acrylic stand for Excel Electronic Pipettor, holds 3 pipettors
P3935L	Carousel stand for Excel Electronic Pipettor, holds 6 pipettors
T01-000131-00	Travel charger and pin for US
T01-000134-00	Travel charger and pin for EU
125-161WP0-02UK	Travel charger and pin for UK
T01-000133-00	Travel charger and pin for AU



Excel™ Electronic Pipettors on acrylic stand



Excel Electronic Pipettors on carousel stand