

Pipet Tips

Standard and filtered

Pipet tips are integral components of any pipetting system. Their shape, material properties and purity have an influence on the pipetting results. These pipet tips are manufactured from virgin polypropylene in clean room conditions and are certified free from bio-contaminants. The range includes standard, sterile, and filtered versions in all standard packaging formats.

Bulk Packs

Pipet tips available in bulks – resealable plastic bags – keeping them safe from contamination. One bag contains 200, 250 or 1,000 tips depending on the tip volume.



Racks

Racked pipet tips delivered in a durable polypropylene box for good stability on the lab bench. One rack contains 96 or 100 tips depending on the tip volume.



Refill system

One refill contains 10 trays with 96 tips each plus an empty rack box. Available only for 10 µL and 200 µL tips with graduation marks.



- ▶ Perfect for use with HTL pipettors for optimal accuracy and precision
- ▶ Compatible with many pipettor brands
- ▶ Strict quality standards (pursuant to ISO 9001) followed from R&D to production and delivery
- ▶ Premium quality virgin polypropylene used as raw material
- ▶ Suitable for sterilization with autoclave (121°C, 15 minutes, w/o dry cycle)

Standard Tips

Cat. No.	Volume (µL)	Pipettor volume range (µL)	Sterile	Packaging	Qty.
301-03-052	0.1 - 10 (short)	0.1 - 2 0.5 - 10	No	Bulk	1,000 tips/bag
301-03-102			No	Rack	96 tips/box, 10 boxes
301-03-152			Yes	Rack	96 tips/box, 10 boxes
301-79-102			No	Refill	96 tips/tray, 10 removable trays*
301-22-052	0.1 - 10 (extended)	0.1 - 2 0.5 - 10	No	Bulk	1,000 tips/bag
301-22-102			No	Rack	96 tips/box, 10 boxes
301-22-152			Yes	Rack	96 tips/box, 10 boxes
301-02-052	1 - 200 (without graduation marks)	2 - 20 5 - 50 10 - 100 20 - 200	No	Bulk	1,000 tips/bag
301-02-152			No	Rack	96 tips/box, 10 boxes
301-02-202			Yes	Rack	96 tips/box, 10 boxes
301-06-052	1 - 200 (graduation marks at 10, 50 and 100 µL)	2 - 20 5 - 50 10 - 100 20 - 200	No	Bulk	1,000 tips / bag
301-06-352			Yes	Rack	100 tips/box, 10 boxes
301-79-052			No	Refill	10 removable trays of 96 tips*
301-16-052	1 - 300 (fine-point)	50 - 250 50 - 300	No	Bulk	1,000 tips/bag
301-16-102			No	Rack	96 tips/box, 10 boxes
301-01-202	100 - 1,000 (bevelled)	100 - 1,000	No	Bulk	1,000 tips/bag
301-01-302			No	Rack	100 tips/box, 10 boxes
301-01-352			Yes	Rack	100 tips/box, 10 boxes
301-71-052	500 - 5,000	500 - 5,000	No	Bulk	250 tips/bag
301-75-052	1,000 - 10,000	1,000 - 10,000	No	Bulk	200 tips/bag

*Refill contains 10 trays with tips plus 2 racks.

Filtered Tips, Rack

Cat. No.	Volume (µL)	Pipettor volume range (µL)	Sterile	Qty.
302-05-102	0.1 - 10 (short)	0.1 - 2 0.5 - 10	Yes	96 tips/box, 10 boxes
302-26-152	0.1 - 10 (extended)	0.1 - 2 0.5 - 10	Yes	96 tips/box, 10 boxes
302-03-102	0.1 - 20 (graduation marks)	2 - 20	Yes	96 tips/box, 10 boxes
302-08-102	1 - 100 (graduation marks)	10 - 100	Yes	96 tips/box, 10 boxes
302-04-102	1 - 200 (graduation marks)	20 - 200	Yes	96 tips/box, 10 boxes
302-12-102	1 - 300 (fine-point)	50 - 250 50 - 300	Yes	96 tips/box, 10 boxes
302-01-102	100 - 1,000 (bevelled)	100 - 1,000	Yes	100 tips/box, 10 boxes

