

for Labnet Pipettors



NOTE: Labnet pipettors are calibrated using non-filtered tips, thus pipettor recalibration might be required when using filter tips.

NOTE: Labnet pipettors are calibrated using non-filtered tips, thus pipettor recalibration might be required when using filter tips.				Labnet BioPette Plus and BioPette A																Labnet BioPette Pro								Labnet Labpette FX																Labnet Excel													
				Channels	Single-channel								Multi-channel				Single-channel				Multi-channel				Single-channel																Single-channel				Multi-channel												
					0.1-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20-200 µL	100-1000 µL	1000-5000 µL	1000-10000 µL	0.5-10 µL	5-50 µL	20-200 µL	50-300 µL	0.1-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20-200 µL	100-1000 µL	0.5-10 µL	5-50 µL	20-200 µL	50-300 µL	2 µL	5 µL	10 µL	20 µL	40 µL	50 µL	70 µL	75 µL	100 µL	120 µL	150 µL	200 µL	250 µL	300 µL	400 µL	500 µL	1,000 µL	0.5-10 µL	2-20 µL	10-200 µL	100-1200 µL	0.5-10 µL	2-20 µL	10-200 µL	100-1200 µL						
	Pipet Tip	Cat. No.	Volume Range																																																						
Labnet	BioFree™	P1000-10	0.5-10 µL	•	•						•				•					•																			•	•																	
		P1000-200	1-200 µL			•	•	•				•	•				•	•	•		•	•			•	•	•	•	•	•	•	•	•	•				•																			
		P1000-300	1-300 µL											•																																											
		P1000-1000	100-1000 µL							•										•																	•	•	•	•	•																
		P3998	1-5 mL							•																																															
		K943	1-10 mL																																																						
Corning®	DeckWorks	4125	0.1-10 µL		•						•				•					•																																					
		4135	0.1-10 µL		•						•				•					•																																					
		4136	1-20 µL			•									•																																										
		4137	1-100 µL				•								•																																										
		4126	1-200 µL			•	•	•				•	•		•	•	•		•	•		•	•		•	•	•	•	•	•	•	•	•	•	•	•																					
		4138	1-200 µL					•				•				•																																									
		4139	1-300 µL					•				•				•																																									
		4128	1-300 µL					•				•	•			•								•																																	
		4129	100-1000 µL						•										•																		•	•	•	•	•																
		4140	100-1000 µL						•										•																																						
	Isotip	4801	0.1-2 µL		•									•						•																																					
		4807	0.2-10 µL		•	•					•			•	•					•																																					
		4808	0.5-10 µL			•					•			•						•																																					
		4821	1-30 µL				•							•						•																																					
4823		1-200 µL					•				•				•								•																																		
4809		100-1000 µL						•										•																																							
Axygen®	MultiRack	MR-10XT*	0.1-10 µL		•					•				•					•																																						
		MRF-10XT*	0.1-10 µL		•					•				•					•																																						
		MRF-20	1-20 µL			•									•																																										
		MRF-100	1-100 µL				•								•																																										
		MR-200	1-200 µL				•	•				•			•	•				•				•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•						
		MRF- 200XT*	1-200 µL				•	•				•			•	•				•				•																																	
		MRF-200NX	1-200 µL				•					•				•				•				•																																	
		MRF- 300XT*	1-300 µL				•					•	•			•				•				•	•																																
		MR-300XT*	1-300 µL				•					•	•			•				•				•	•																																
		MR-300	1-300 µL				•					•	•			•				•				•	•																																
		MR-1000XT*	50 - 1250 µL					•											•																																						
		MRF-1000XT*	50 - 1250 µL					•											•																																						
	Universal Fit	T-300	0.5-10 µL		•	•					•			•	•					•																																					
		TF-300	0.5-10 µL		•	•					•			•	•					•																																					
		TXL-10	0.5-20 µL			•					•				•					•																																					
		TF-20	0.5-20 µL				•								•																																										
		TF-50	0.5-50 µL								•									•																																					
		TF-100	1-100 µL				•								•																																										
		TR-222	1-200 µL				•	•				•				•	•			•	•			•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•				
		TF-200	1-200 µL					•				•								•				•																																	
		T-350	1-300 µL								•	•								•	•			•	•																																
		TF-350	1-300 µL					•				•	•							•	•			•	•																																
		T-1000	100-1000 µL						•											•																																					
		TF-1000	100-1000 µL						•											•																																					
TF-1200	100-1100 µL						•											•																																							
T-5000	1-5 mL							•																																																	
TF-10ML	1-10 mL								•																																																

*Extended length tip

Corning Incorporated
Life Sciences
www.corning.com/lifesciences

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India
t 91 124 4604000

Japan
t 81 3-3586 1996

Korea
t 82 2-796-9500

Singapore
t 65 6572-9740

Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com

LATIN AMERICA
grupolA@corning.com

Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2023, 2024 Corning Incorporated. All rights reserved. 1/24 CLSLN-EQ-023 REV1

