

Labnet FastPette™ Pro Pipet Controller

Lightweight, ergonomically shaped, designed for work convenience

The FastPette Pro is a lightweight and ergonomically designed automatic pipet controller ideal for liquid handling with 0.5 to 100 mL glass or plastic serological pipets. Conveniently positioned switches allow toggling of operation modes and aspiration speeds for high precision handling to suit different pipet volumes and liquid viscosities. Aspirating and dispensing speed is controlled by two comfortable finger triggers. Liquids can be dispensed by gravity or by blow out using the motorized pump.

To protect the device against overfilling, the FastPette Pro is equipped with a 0.2 µm PTFE filter. The PTFE filter, pipet holder, and nose piece are fully autoclavable and can be individually replaced as needed.

Features

- ▶ Triggers convenient for one or two finger control of speed and mode
- ▶ LCD display for convenient visualization of speed, modes, and battery status
- ▶ Easy access to the battery compartment. Li-ion battery allows cordless work for up to 6 hours
- ▶ 1 year warranty

Ordering Information

Cat. No.	Description
P2004	FastPette Pro Pipet Controller, with 0.2 µm filter, Li-ion battery, universal charger with set of adapters, and wall hanger

Accessories

Cat. No.	Description	Qty/Cs
SP9143	PTFE autoclavable filter, 0.2 µm	5
SP9144	PTFE autoclavable filter, 0.45 µm	5
SP29052	Nose piece	1
SP29054	Silicone pipet holder	1
SP29056	Charging stand, 2-positions	1
SP9029	Wall hanger	1
SP29096	Battery set, Li-ion	1
SP29091	Battery cover	1
SP29100	Universal power supply, 9V with set of adapters (EU, USA, GBR, KOR, and AUS)	1

