

	Finntip (sterile and non-sterile)															
	10 Micro	10 Flex micro	20 Micro	50 Micro	250 Univ	200 Ext	200 Flex	300	300 Flex	1000 Ext	1000	1000 Flex	1200 Flex	5 ml	10ml	10ml Flex Ext
FP F2 0.2 - 2 µl	•	•	•	•												
FP F2 0.5 - 5 µl	•	•	•	•												
FP F2 1 - 10 µl micro	•	•	•	•												
FP F2 1 - 10 µl					•	•	•	•	•							
FP F2 2 - 20 µl micro			•	•												
FP F2 2 - 20 µl					•	•	•	•	•							
FP F2 5 - 50 µl micro				•												
FP F2 5 - 50 µl					•	•	•	•	•							
FP F2 10 - 100 µl					•	•	•	•	•							
FP F2 20 - 200 µl					•	•	•	•	•							
FP F2 100 - 1000 µl										•	•	•	•			
FP F2 0.5 - 5 ml														•		
FP F2 1 - 10 ml															•	•
FP F2 Fixed Volume 1 µl	•	•	•	•												
FP F2 Fixed Volume 5 µl	•	•	•	•												
FP F2 Fixed Volume 10 µl					•	•	•	•	•							
FP F2 Fixed Volume 20 µl					•	•	•	•	•							
FP F2 Fixed Volume 25 µl					•	•	•	•	•							
FP F2 Fixed Volume 50 µl					•	•	•	•	•							
FP F2 Fixed Volume 100 µl					•	•	•	•	•							
FP F2 Fixed Volume 200 µl					•	•	•	•	•							
FP F2 Fixed Volume 250 µl										•	•	•	•			
FP F2 Fixed Volume 500 µl										•	•	•	•			
FP F2 Fixed Volume 1000 µl										•	•	•	•			
FP F2 Fixed Volume 2000 µl														•		
FP F2 Fixed Volume 3000 µl														•		
FP F2 Fixed Volume 5000 µl														•		
FP F2 Fixed Volume 10000 µl															•	•
FP F2 8-Ch 1-10 µl	•	•	•	•												
FP F2 12-Ch 1-10 µl	•	•	•	•												
FP F2 8-Ch 5-50 µl					•	•	•	•	•							
FP F2 12-Ch 5-50 µl					•	•	•	•	•							
FP F2 8-Ch 10-100 µl					•	•	•	•	•							
FP F2 12-Ch 10-100 µl					•	•	•	•	•							
FP F2 8-Ch 30-300 µl								•	•							
FP F2 12-Ch 30-300 µl								•	•							
FP F2 16-Ch 1-10 µl			•	•												
FP F2 16-Ch 5-50 µl				•												

	Finntip Filter (sterile)																						
	10 Micro	10 Flex micro	20 Micro	50 Micro	10 Univ	20 Univ	30 Univ	30 Flex	100 Ext	100 Univ	100 Flex	200 Ext	200 Flex	200 Univ	300	300 Flex	1000 Ext	1000	1000 Flex	1200 Flex	5 ml	10ml	10ml Flex Ext
FP F2 0.2 - 2 µl	•	•	•	•																			
FP F2 0.5 - 5 µl	•	•	•	•																			
FP F2 1 - 10 µl micro	•	•	•	•																			
FP F2 1 - 10 µl					•	•	•	•															
FP F2 2 - 20 µl micro			•	•																			
FP F2 2 - 20 µl						•	•	•															
FP F2 5 - 50 µl micro				•																			
FP F2 5 - 50 µl									•	•	•	•	•	•	•	•							
FP F2 10 - 100 µl									•	•	•	•	•	•	•	•							
FP F2 20 - 200 µl												•	•	•	•	•							
FP F2 100 - 1000 µl																	•	•	•	•			
FP F2 0.5 - 5 ml																					•		
FP F2 1 - 10 ml																						•	•
FP F2 Fixed Volume 1 µl	•	•	•	•																			
FP F2 Fixed Volume 5 µl	•	•	•	•																			
FP F2 Fixed Volume 10 µl					•	•	•	•															
FP F2 Fixed Volume 20 µl						•	•	•															
FP F2 Fixed Volume 25 µl							•	•															
FP F2 Fixed Volume 50 µl									•	•	•	•	•	•	•	•							
FP F2 Fixed Volume 100 µl									•	•	•	•	•	•	•	•							
FP F2 Fixed Volume 200 µl												•	•	•	•	•							
FP F2 Fixed Volume 250 µl																	•	•	•	•			
FP F2 Fixed Volume 500 µl																	•	•	•	•			
FP F2 Fixed Volume 1000 µl																	•	•	•	•			
FP F2 Fixed Volume 2000 µl																						•	
FP F2 Fixed Volume 3000 µl																						•	
FP F2 Fixed Volume 5000 µl																						•	
FP F2 Fixed Volume 10000 µl																							•
FP F2 8-Ch 1-10 µl	•	•	•	•																			
FP F2 12-Ch 1-10 µl	•	•	•	•																			
FP F2 8-Ch 5-50 µl									•	•	•	•	•	•	•	•							
FP F2 12-Ch 5-50 µl									•	•	•	•	•	•	•	•							
FP F2 8-Ch 10-100 µl									•	•	•	•	•	•	•	•							
FP F2 12-Ch 10-100 µl									•	•	•	•	•	•	•	•							
FP F2 8-Ch 30-300 µl																•	•						
FP F2 12-Ch 30-300 µl																•	•						
FP F2 16-Ch 1-10 µl			•	•																			
FP F2 16-Ch 5-50 µl				•																			

Tip Ordering Information
Bestellinformation: Pipettenspitzen
Renseignements pour commander des cônes
Información para pedidos del portapuntas
チップの注文について

Code	Finntip	Volume	Qty
* 9400300	10 Micro	0,2-10 µl	10x96/rack
9400303	10 Micro sterile	0,2-10 µl	10x96/rack
* 94060100	Flex 10	0,2-10 µl	10x96/rack
94060103	Flex 10 sterile	0,2-10 µl	10x96/rack
* 9400610	20 Micro	0,2-20 µl	10x384/rack
9400613	20 Micro sterile	0,2-20 µl	10x384/rack
* 9400370	50 Micro	0,2-50 µl	10x384/rack
9400373	50 Micro sterile	0,2-50 µl	10x384/rack
* 9400130	200 Ext	5-200 µl	10x96/rack
9400133	200 Ext sterile	5-200 µl	10x96/rack
* 94060310	Flex 200	1-200 µl	10x96/rack
94060313	Flex 200 sterile	1-200 µl	10x96/rack
* 9400260	250 Univ.	0,5-250 µl	10x96/rack
9400263	250 Univ. sterile	0,5-250 µl	10x96/rack
* 9401250	300	5-300 µl	10x96/rack
9401253	300 sterile	5-300 µl	10x96/rack
* 94060510	Flex 300	5-300 µl	10x96/rack
94060513	Flex 300 sterile	5-300 µl	10x96/rack
* 9401110	1000	100-1000 µl	10x96/rack
9401113	1000 sterile	100-1000 µl	10x96/rack
* 94060710	Flex 1000	100-1000 µl	10x96/rack
94060713	Flex 1000 sterile	100-1000 µl	10x96/rack
* 9402070	5 ml	1-5 ml	5x54/rack
9402073	5 ml sterile	1-5 ml	5x54/rack
* 9402180	10 ml	1-10 ml	5x24/rack
9402183	10 ml sterile	1-10 ml	5x24/rack
94060970	Flex 10ml Ext	1-10 ml	100 pcs/bag
94060973	Flex 10ml Ext sterile	1-10 ml	50 pcs/bag

*** Also available in bags**

Code	Finntip Filter	Volume	Qty
94052000	10 Micro sterile	0,2-10 µl	10x96/rack
94056980	Flex 10 sterile	0,2-10 µl	10x96/rack
94052100	10 sterile	0,5-10 µl	10x96/rack
94052020	20 Micro sterile	0,2-20 µl	10x384/rack
94052150	20 sterile	0,5-20 µl	10x96/rack
94052160	30 sterile	0,5-30 µl	10x96/rack
94056510	Flex 30 sterile	5-30 µl	10x96/rack
94052060	50 Micro sterile	0,2-50 µl	10x384/rack
94052200	100 µl sterile	0,5-100 µl	10x96/rack
94052310	100 µl Ext sterile	5-100 µl	10x96/rack
94056520	Flex 100 sterile	5-100 µl	10x96/rack
94052300	200 µl sterile	0,5-200 µl	10x96/rack
94052320	200 µl Ext sterile	5-200 µl	10x96/rack
94056380	Flex 200 sterile	1-200 µl	10x96/rack
94052350	300 µl sterile	5-300 µl	10x96/rack
94056580	Flex 300 sterile	5-300 µl	10x96/rack
94052410	1000 µl sterile	100-1000 µl	10x96/rack
94052430	1000 µl Ext sterile	100-1000 µl	5x96/rack
94056710	Flex 1000 sterile	100-1000 µl	10x96/rack
94052550	5 ml sterile	1-5 ml	5x54/rack
94052800	10 ml sterile	1-10 ml	5x24/rack
94056970	Flex 10 ml Ext sterile	1-10 ml	50/bag