Eppendorf Pipette Holder System

Well within reach: the Eppendorf Pipette Holder System helps you to organize your lab
Save valuable bench space with the Eppendorf Pipette Holder System! With space for up to six single- or multi-channel pipettes, you can easily store all of your Eppendorf pipettes and multi-dispensers within easy reach. Our flexible carousels, stands and wall mount devices fit current Eppendorf manual liquid handling instruments as well as most of their predecessors.

New



> Pipette Carousel 2: Rotatable holder for six mechanical pipettes
> Charger Carousel 2: Carousel with intelligent charging electronics, capacity for up to six instruments (charges four instruments at a time and all six in turn) and magnetic power connector for fast and easy connection



> With Charger Stand 2, the Eppendorf Xplorer®/Xplorer plus or Multipette® E3/E3x is always fully charged and ready for use. Pipette Stand 2 carries one pipette or multi-dispenser
> Large rubber feet protect stands and carousels from liquids spilled on lab bench
> Charger Stand 2 with cord wrap functionality for more convenience

New



DESIGN AWARD*

reddot award 2018 winner

* Golden Pin Design Award Winner 2017

> Place your pipette and Multipette where you need them: Mounted to a wall, on the shelf above your bench or inside a biological safety cabinet
> Pictograms on the reversible holders allow for easy assignment of pipettes to their respective holder

Ordering information

Description	Order no.
New **Pipette Carousel 2, white**, with 6 holders for Eppendorf Research® 3 pipettes, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116000236
Pipette Carousel 2, with 6 holders for Eppendorf Research®, Research plus, Reference, Reference 2, Biomaster® pipettes or epMotion® Tip Tools, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116000015
Charger Carousel 2, with 6 charger shells for Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipettes and mains/power adapter, additional charger shells and holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116000023
Charger Stand 2, with charger shell for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, operated with mains/power adapter supplied with Eppendorf electronic pipettes and dispensers	3116000031
Charger Stand 2, with charger shell for one Multipette® E3/E3x or Multipette® stream/Xstream electronic multi-dispenser, operated with mains/power adapter supplied with Eppendorf electronic pipettes and dispensers	3116000040
Pipette Stand 2, with pipette holder for one Multipette® M4 multi-dispenser, without charging function, additional pipette holders, compatible with other Eppendorf pipettes and dispensers, are sold separately	3116000058
New **Pipette Holder 2, white**, for one Eppendorf Research® 3 pipette, for Pipette Carousel 2, Charger Carousel 2, Pipette Stand 2 or wall mounting, sticky tape included	3116000295
Pipette Holder 2, for one Multipette® E3/E3x or Multipette® stream/Xstream electronic multi-dispenser, for Pipette Carousel 2, Charger Carousel 2, Pipette Stand 2 or wall mounting, sticky tape included, without charging function	3116000139
Pipette Holder 2, for one Eppendorf Research®, Research plus, Eppendorf Reference®, Reference 2, Biomaster® pipette or epMotion® Tip Tool, for Pipette Carousel 2, Charger Carousel 2, Pipette Stand 2 or wall mounting, sticky tape included, without charging function	3116000112
Pipette Holder 2, for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, for Pipette Carousel 2, Charger Carousel 2, Pipette Stand 2 or wall mounting, sticky tape included, without charging function	3116000120
Pipette Holder 2, for one Multipette® M4 multi-dispenser, for Pipette Carousel 2, Charger Carousel 2, Pipette Stand 2 or wall mounting, sticky tape included, without charging function	3116000147
Charger Shell 2, for one Eppendorf Xplorer®, Xplorer plus or Xplorer plus Move It® electronic pipette, for Charger Carousel 2 or Charger Stand 2, with charging function	3116602007
Charger Shell 2, for one Multipette® E3/E3x or Multipette® stream/Xstream electronic multi-dispenser, for Charger Carousel 2 or Charger Stand 2, with charging function	3116603003